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BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

Do you need more operating manuals and translations?
Please refer to <http://www.brand.de/om> or use the following
Quick Response Code:

The original operating manual is written in German. Other languages are translations of the original operating manual.

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The manufacturer is solely responsible for issuing this declaration of conformity with regard to meeting the essential requirements and preparing the technical documentation.

Hersteller / Manufacturer	Adresse / Address
BRAND GMBH + CO KG	Otto-Schott-Str. 25 97877 Wertheim Germany

Wertheim, 25. Januar 2023 / 25 January 2023

Patrick Ziemeck
Managing Director Technology

i.A. Siegfried Ott
Regulatory Compliance

UK Declaration of Conformity

Product name: Pipette leak tester PLT unit

Product variations: -

Cat.-No. 703970

The product described above meets the basic requirements that are specified in the harmonization legislation listed below:	Applied harmonized standards:
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 (SI 2012 No. 3032)	EN IEC 63000:2018
Electromagnetic Compatibility Regulations 2016 (SI 2016 No. 1091, as amended by SI 2019 No. 696)	EN 61326-1:2013
Other applied standards:	Scope:
EN 61010-1:2010, EN 61010-1:2010/A1:2019/AC:2019-04, EN 61010-1:2010/A1:2019 EN/IEC 62368-1:2014 (applied to third-party power supply)	Electrical and functional safety Electrical and functional safety

The manufacturer is solely responsible for issuing this declaration of conformity with regard to meeting the essential requirements and preparing the technical documentation.

Manufacturer	Address
BRAND GMBH + CO KG	Otto-Schott-Str. 25 97877 Wertheim Germany

Wertheim, 23 August 2022

Patrick Ziemeck
Managing Director Technology

i.A. Siegfried Ott
Regulatory Compliance

This document declares the accordance with the named harmonized regulations, but does not assure specific properties.

Declaration of conformity – China RoHS 2

BRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in BRAND products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB/T 26572 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

+ Lead (Pb):	0.1%	+ Hexavalent chromium (Cr(+VI)):	0.1%
+ Mercury (Hg):	0.1%	+ Polybrominated biphenyls (PBB):	0.1%
+ Cadmium (Cd):	0.01%	+ Polybrominated diphenyl ether (PBDE):	0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for BRAND instruments is 40 years.

Material Content Declaration for BRAND Products

部件名称 Part name	有毒有害物质或元素 Hazardous substances					
	铅 Pb	汞 Hg	镉 Cd	六价铬 Cr(+VI)	多溴联苯 PBB	多溴二苯醚 PBDE
包装 / Packaging	0	0	0	0	0	0
塑料外壳 / 组件 Plastic housing / parts	0	0	0	0	0	0
电池 / Battery	0	0	0	0	0	0
玻璃 / Glass	0	0	0	0	0	0
电子电气组件 Electrical and electronic parts	X	0	0	0	0	0
金属外壳 / 组件 Metal housing / parts	X	0	0	0	0	0
电机 / Motor	X	0	0	0	0	0
配件 / Accessories	X	0	0	0	0	0

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This table is created according to SJ/T 11364-2014.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Note: Batteries, glassware and accessories might not be content of the enclosed device and/or may have their own EFUP-marking and/or might be maintaining parts with changing EFUP-marking.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by BRAND may enter into further devices or can be used together with other appliances. With these third party products and appliances in particular, please note the EFUP labeled on these products. BRAND will not take responsibility for the EFUP of those products and appliances.

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Table of contents

1 Introduction	49
1.1 Package contents.....	49
1.2 Terms of use	49
2 Safety instructions.....	50
2.1 Safety regulations.....	50
2.2 Function and operating limits.....	50
2.3 Excluded uses.....	51
2.4 Operating limits	51
2.5 Technical data.....	51
3 Operating elements.....	52
4 The first steps	53
4.1 Powering on the device	53
4.2 Self-test	53
5 Performing the test	58
5.1 Test preparation	58
5.2 Changing the adapter and filter	60
5.3 Testing single-channel pipettes.....	62
5.4 Testing multichannel pipettes	65
6 Switching the device off	68
7 Make settings	69
7.1 Select the language	69
7.2 Set the contrast.....	70
7.3 Inspection date	71
7.4 Auto Power Off Time.....	72
7.5 Pressure units (hPa, mbar, torr)	73
7.6 Factory settings.....	74
8 Firmware update	76
9 Correlation table.....	77
10 Troubleshooting	78
11 Product markings	80
12 Ordering information	81
13 Calibration service	84
14 Warranty	85
15 Repairs.....	86
15.1 Sending for repair	86
16 Disposal	88

1 Introduction

1.1 Package contents

The package contains the PLT unit, one adapter each for testing air-displacement pipettes with tips (mounted) and without tips, two caps, three replacement PE filters for the pipette adapters, an AC adapter, a quality certificate, and this Operating Manual.

1.2 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-do-date versions of the operating manual on our website: www.brand.de.

1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2.2 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

2 Safety instructions

2.1 Safety regulations

Please read carefully!

The laboratory device can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- All users must read this operating manual and heed the information contained in it before using the device.
- Follow the general safety rules and regulations; e.g., wear protective clothing, eye protection, and safety gloves.
- Follow standard safe laboratory practices and procedures when working with infectious or hazardous samples.
- Do not operate the device in a potentially explosive atmosphere.
- Only use the device for leak-testing pipettes that operate according to the air displacement principle. Ensure that the device is not used for any of the excluded uses described in Operating exclusions! When in doubt, contact the manufacturer or vendor.
- Use only in a manner that poses no risk to the user or other people.
- When testing air-displacement pipettes with mounted tips, employ only unused tips.
- Never use force.
- Only use original replacement parts. Do not make any modifications to the device.
- Repairs and work inside the device should only be conducted by authorized electricians.
- Always ensure that the device is in good condition before using it. If there are indications of a problem with the device, interrupt the test immediately and follow the instructions in Troubleshooting. If necessary, contact the manufacturer.
- The AC adapter should be protected from moisture and may only be used with this device.

WARNING

Improper handling of the device (short circuit, mechanical damage, overheating, wrong AC adapter, etc.) can, in extreme cases, lead to fire or explosion of the device.

2.2 Function and operating limits

The BRAND leak testing unit is a differential pressure measuring device that determines leak rates in air-displacement pipettes. Static and dynamic testing permit the localization of leaks. Dynamic testing is preferable.

Testing is possible both with tips (for testing the instrument/tip interface and the entire system) and without tips (for testing the instrument).

The leak rate is directly related to the volume loss in the pipette.

The leak rate Q_L is a measure of the amount of substance (mass) flowing through a leak per unit of time.

The leak rate is the ratio of the pV value of a gas to the period of time during which the gas flows through the leak. The pV value is the product of the pressure and the volume of a specific amount of gas at a given temperature. For an ideal gas, the pV value at a given temperature is a measure of the amount of substance or the mass of the gas. The leak rate is dependent on the type of gas, the absolute pressure, the pressure difference, and the temperature.

In the pipette tests, suitable units for the leak rate Q_L are hPa mL/s. A leak rate of 1 hPa mL/s at an air pressure of 1000 hPa corresponds to a volume loss of 1 μ L/s.

NOTICE

The leak testing unit cannot be used in lieu of regular gravimetric testing within the scope of measuring instrument monitoring activities. The device should be used in daily checks of the pipettes as a safeguard for the periods between calibrations.

2.3 Excluded uses

- The user is responsible for checking whether the device is suitable for a given application.
- Do not use the device to aspirate liquids! Avoid aggressive vapors (corrosion hazard)! Further information is available on request.
- Do not evacuate hollow vessels that are not specified for vacuum use; otherwise, the vessel may burst!

2.4 Operating limits

The device offers a rapid and simple test for leaks in air-displacement pipettes within the following limits:

- Operating temperature 15 °C to 35 °C
- Relative humidity max. 80%
- Elevation at place of use max. 2000 meters above sea level
- Storage temperature 0 °C to 50 °C

NOTICE


If the device is relocated, variations in temperature and humidity during transport and in the new environment can affect the measurements.

To avoid error messages, allow the PLT unit to acclimate in the testing room for a period of 2 hours before using it!

2.5 Technical data

PLT unit	Input: DC 6.5 V, 800 mA
AC adapter	Input: AC 100 V ... 240 V, 50/60 Hz, 0,5 A Output: DC 6.5 V, 800 mA

3 Operating elements

- 1 Single-channel adapter
- 2 Display
- 3 "Start/Execute" button
- 4 Select knob
- 5 LED indicator
- 6 "Back" button
- 7 AC adapter socket
- 8 USB port

"Back" button




The function of the "Back" button differs depending on whether it is pressed for a short or a long time. Pressing the button for a short time leads back from a submenu item to the next higher level of the menu hierarchy. Pressing the button for a long time leads back directly to the "Ready" start menu.

Select knob

In many cases, you can press the left "Start/Execute" button or the central knob to confirm a selection. For the sake of consistency, the "Start/Execute" button on the left is used for confirmation in this operating manual.

4 The first steps

4.1 Powering on the device

What do I need to do?	How is that going?	Which button?	What does the display show?
<p>1. Connect the power supply</p> <p>Use only original power supply unit! The power supply unit and power supply plug on the device must be freely accessible and must be able to be disconnected from the mains at any time.</p>	<p>Plug the power supply into the socket on the back of the device, without pressing any button on the device.</p> <p>When the power supply is connected, the display shows the software version number, the date of creation and the BRAND logo.</p> <p>After 3 s, the display jumps to the 'Language' display during the first start-up.</p>		 <p>Firmware V1.40.011 3.1.2022</p> <p>BRAND</p>
2. Select a language	Select the language by turning the selection button.		 <p>Language Deutsch English Español o.k.</p>
3. Confirm the language	Press the left confirmation button. The 'Settings' menu appears in the selected language. The Back button takes you to the main menu.		 <p>Settings Language Contrast Inspection date Auto-Power-Off Time select back</p>

NOTICE


After commissioning, you can either perform a self-test or start the pipette test directly.

4.2 Self-test

4.2.1 Quick test

The self-test determines if there are any leaks in the BRAND PLT unit. Two different tests are available. The Quick Test lasts only about 2 minutes and should be performed on a daily basis. The Full Test

lasts about 12 minutes. It includes the additional differentiated testing of the internal tanks and should be performed on a weekly basis.


What to do	How to do it	Button to use	Display
1. Call up the Main Menu	Press or turn the select knob to call up the Main Menu.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <div style="border: 1px solid white; padding: 2px; margin-bottom: 2px;">Pipette type</div> <p>Self-test Device OFF Settings</p> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> Selection back </div> </div>
2. Set up the self-test	Turn the select knob to the "Self-Test" setting.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <p>Pipette type</p> <div style="border: 1px solid white; padding: 2px; margin-bottom: 2px;">Self-test</div> <p>Device OFF Settings</p> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> Selection back </div> </div>
3. Confirming the self-test	Press the "Execute" button. The "Self-Test" submenu is called up. The last test performed is shown with a frame around it.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <p>Quick Test</p> <div style="border: 1px solid white; padding: 2px; margin-bottom: 2px;">Full test</div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> Selection back </div> </div>
4. Select the Quick Test	Turn the select knob to the "Quick Test" setting.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Self Test</p> <div style="border: 1px solid white; padding: 2px; margin-bottom: 2px;">Quick test</div> <p>Full test</p> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div> </div>
5. Confirm the Quick Test	Press the "Start/Execute" button. The start menu is called up.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Self-test</p> <p>Please insert the blanking plug</p> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> Selection back </div> </div>
6. Seal the system	Press the supplied cap firmly into the adapter (for the multichannel adapter, seal all four channels).		
7. Start the Quick Test	Press the "Start/Execute" button. The five components – sensor, pump, valve, filter, and measurement chamber – are tested in succession.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Self-test</p> <div style="border: 1px solid white; padding: 2px; margin-bottom: 2px;"> <div style="width: 100%; height: 10px; background-color: white;"></div> </div> <p style="text-align: center;">-0.10 hPa 1:sensor</p> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> back </div> </div>


What to do	How to do it	Button to use	Display
	The progress bar on the display indicates the test status. The LED indicator alternately blinks red and green during the test.		
8. Aerate the system	Once the component test has been completed, the display shows a prompt to remove the cap.		<div style="background-color: #333; color: white; padding: 5px; text-align: center;">Self-test</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Please pull the dummy plug</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Selection back</div>
	Remove the cap and press the "Start/Execute" button. The system is aerated; the display shows the current pressure in the measurement chamber.		
9. End the self-test	Once the system has been aerated, the display automatically shows the result of the test. If an error occurs, the system stops the procedure, an error code is displayed (see Troubleshooting, p. 78), and the system is aerated. To repeat the test, press the "Start/Execute" button. To go back to the Main Menu, press the "Back" button.		<div style="background-color: #333; color: white; padding: 5px; text-align: center;">Self-test</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Passed</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Restart back</div>
			<div style="background-color: #333; color: white; padding: 5px; text-align: center;">Self-test</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Error 1XXX see manual</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Restart back</div>

4.2.2 Full test

The self-test determines if there are any leaks in the BRAND PLT unit. Two different tests are available. The Quick Test lasts only about 2 minutes and should be performed on a daily basis. The Full Test lasts about 12 minutes. It includes the additional differentiated testing of the internal tanks and should be performed on a weekly basis.

What to do	How to do it	Button to use	Display
1. Call up the Main Menu	Press or turn the select knob to call up the Main Menu.		<div style="background-color: #333; color: white; padding: 5px; text-align: center;">Main menu</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Pipette type</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Self-test</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Device OFF</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Settings</div> <div style="background-color: #333; color: white; padding: 5px; text-align: center;">Selection back</div>

What to do	How to do it	Button to use	Display
2. Set up the self-test	Turn the select knob to the "Self-Test" setting.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <p>Pipette type</p> <p style="border: 1px solid white; padding: 2px;">Self-test</p> <p>Device OFF</p> <p>Settings</p> </div> <div style="display: flex; justify-content: space-between; border-top: 1px solid #ccc; margin-top: 5px;"> Selection back </div>
3. Confirming the self-test	Press the "Execute" button. The "Self-Test" submenu is called up. The last test performed is shown with a frame around it.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <p>Quick Test</p> <p style="border: 1px solid white; padding: 2px;">Full test</p> </div> <div style="display: flex; justify-content: space-between; border-top: 1px solid #ccc; margin-top: 5px;"> Selection back </div>
4. Select the Quick Test	Turn the select knob to the "Full Test" setting.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <p>Quick Test</p> <p style="border: 1px solid white; padding: 2px;">Full test</p> </div> <div style="display: flex; justify-content: space-between; border-top: 1px solid #ccc; margin-top: 5px;"> Selection back </div>
5. Confirm the Quick Test	Press the "Start/Execute" button. The start menu is called up.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Self-test</p> <p>Please insert the blanking plug</p> </div> <div style="display: flex; justify-content: space-between; border-top: 1px solid #ccc; margin-top: 5px;"> Selection back </div>
6. Seal the system	Press the supplied cap firmly into the adapter (for the multichannel adapter, seal all four channels).		
7. Start the Quick Test	Press the "Start/Execute" button. The five components – sensor, pump, valve, filter, and measurement chamber – are tested in succession. The progress bar on the display indicates the test status. The LED indicator alternately blinks red and green during the test.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Self-test</p> <div style="border: 1px solid white; width: 100%; height: 15px; margin: 5px 0;"></div> <p style="text-align: center; margin: 0;">-0.10 hPa 1:sensor</p> </div> <div style="display: flex; justify-content: space-between; border-top: 1px solid #ccc; margin-top: 5px;"> back </div>

What to do	How to do it	Button to use	Display
8. Aerate the system	Once the component test has been completed, the display shows a prompt to remove the cap.		<p>Self-test</p> <p>Please pull the dummy plug</p> <p>Selection back</p>
	Remove the cap and press the "Start/Execute" button. The system is aerated; the display shows the current pressure in the measurement chamber.		<p>Self-test</p> <p>-20.10 hPa</p> <p>back</p>
9. End the self-test	Once the system has been aerated, the display automatically shows the result of the test. If an error occurs, the system stops the procedure, an error code is displayed (see Troubleshooting, p. 78), and the system is aerated. To repeat the test, press the "Start/Execute" button.		<p>Self-test</p> <p>Passed</p> <p>Restart back</p>
			<p>Self-test</p> <p>Error 1XXX see manual</p> <p>Restart back</p>
	To go back to the Main Menu, press the "Back" button.		

5 Performing the test

5.1 Test preparation

Before the test, you must decide whether a dynamic or static test is to be performed.

Dynamic or static test?

Dynamic test

With the **dynamic test**, you can not only detect static leaks but also determine whether a defective (contaminated, scratched) piston is the cause of a leak.

Over the course of the measurement, you should press the plunger button down completely and evenly 2–3 times. The associated piston movement allows piston faults to be recognized.

For 5 mL and 10 mL pipettes, the button should be pressed and released relatively slowly.

In the dynamic testing of electronic pipettes, the piston should be in the uppermost position at the beginning and should be brought back to that position before the end.

Static test

In the **static test**, in contrast, the plunger button is not pressed, i.e., the piston is not moved, during the test.

In electronic pipettes, the piston can be positioned in the uppermost position or at the first

Leaks occurring during the piston stroke thus cannot be identified in this case.

We recommend that the dynamic test be performed.

Test with or without pipette tip?

The test can be performed with or without the pipette tip.

We recommend that the test be performed with a mounted unused tip so that the entire pipette system is tested.

If a leak is identified, you can repeat the test without a tip to determine whether the leak arises from the tip cone/tip coupling region. To perform the test without a tip, you must first change the adapter; see Changing the adapter and filter, p. 60.

Hold the pipette firmly during the test.

Testing of single- or multichannel pipettes

Changing the adapter

A total of three adapters are available:

Adapter 1:

Already premounted; for testing single-channel pipettes up to 10 mL with mounted pipette tips

Adapter 2:

Included; for testing single-channel pipettes without pipette tips

Adapter 3:

Available separately (optional); for testing multichannel pipettes with and without pipette tips, with four channels always being tested simultaneously. The 4-channel adapter can be used for 4-, 8-, and 12-channel pipettes.

If a leak is detected, three channels can be closed with caps. The leaking channel can be identified through subsequent single-channel testing.


During multichannel pipette testing, the tips can be allowed to remain in the adapter when the device is disconnected. They can then simply be removed by hand or left inserted for the next test.

PE Filter

NOTICE

Important!

To protect the device from moisture and contamination, an easily removable PE filter is provided at the bottom of the silicone adapter.


5.2 Changing the adapter and filter

a.


a. Unscrew the adapter lock nut.

b.


b. Remove the pipette adapter.

c.


c. If the PE filter is clogged with moisture or contaminants, it must be changed. In this case, carefully squeeze the PE filter out from the bottom end of the silicone adapter.

d.


d. Insert a new filter.

e.

e. Insert the selected pipette adapter into the adapter housing (Adapter 2 for pipettes without tips shown here).

f.

f. Screw the lock nut in tightly.

g. After changing the adapter, perform the Quick Test.

Changing the 4-channel adapter

The 4-channel adapter consists of an adapter with an integrated lock nut. The entire ensemble is exchanged.

5.3 Testing single-channel pipettes

Example: 100 µL air-displacement pipette, dynamic test

Before performing the test, set the pipette to the nominal volume. Mount a new unused disposable tip, and insert the pipette with the tip into the previously selected and attached Adapter 1. Select an adapter; see Test preparation, p. 58.


What to do	How to do it	Button to use	Display
1. Call up the Main Menu	Press or turn the select knob to call up the Main Menu.		Main menu Pipette type Self-test Device OFF Settings Selection back
2. Set up the self-test	Turn the select knob to the "Pipette Type" setting.		Main menu Pipette type Self-test Device OFF Settings Selection back
3. Confirm the pipette type	Press the "Execute" button. The "Channels" submenu is called up.		Channels Single Channel Multichannel select back
4. Select the Quick Test	Turn the select knob to the "Single Channel" setting.		Channels Single Channel Multichannel select back
5. Confirm the single-channel pipette	Press the "Execute" button. The "Nominal Volume" menu is called up.		Nominal Volume 25 µl 50 µl 100 µl 200 µl select back
6. Select the pipette volume	Turn the select knob to the "100 µl" setting.		Nominal Volume 25 µl 50 µl 100 µl 200 µl select back

What to do	How to do it	Button to use	Display
7. Confirm the pipette volume	Press the "Start/Execute" button. The start submenu is called up.		<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">ready</p> <p style="margin: 0;">SC - 100µl</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
8. Start the test	Press the "Start/Execute" button. The test is begun, and the test pressure is set (preparing...).		<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Test</p> <p style="margin: 0;">preparing...</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
Press the pipette plunger button 2-3 times	Shortly before preparation is finished, the LED light lights up once in red and once in green. Then the progress bar appears and the system automatically starts the measurement. Press the plunger button down evenly 2-3 times during the measurement. The piston must return to the initial position before the measurement is finished! The progress bar indicates the remaining time and the currently achieved vacuum. The LED indicator blinks orange.		<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Test</p> <div style="background-color: white; width: 100%; height: 15px; margin: 5px 0;"></div> <p style="margin: 0;">-121.21 hPa</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div> <p style="margin-top: 10px;">(pressure value shown as an example)</p> <div style="text-align: center; margin-top: 10px;">  <p style="margin: 0;">LED orange</p> </div>
9. End the test	After the test is finished, the display indicates the leak rate that was determined and whether the system passed the test. To repeat the test, press the "Start/Execute" button. To go back to the Main Menu, turn the select knob.	<div style="text-align: center; margin-bottom: 20px;">  <p style="margin: 0;">LED green</p> </div> <div style="text-align: center;">  <p style="margin: 0;">LED red</p> </div>	<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">PASS</p> <p style="margin: 0;">SC – 100 µl Q_L: 0.04 hPa*m/s</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> start </div> <p style="margin-top: 10px;">(leak rate shown as an example)</p> <div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">FAIL</p> <p style="margin: 0;">SC– 100 µl Q_L: 4.6 hPa*m/s</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> start </div> <p style="margin-top: 10px;">(leak rate shown as an example)</p>

NOTICE

You can use the correlation table to estimate the missing volume.

What does the vertical progress bar on the display show?

A mark is added to the progress bar. It indicates the threshold value. It corresponds to approx. $\frac{1}{4}$ of the volume tolerance according to ISO 8655-2.

A: The progress bar is filled in all the way to the bottom: The pipette is leak-tight.

B: The progress bar is filled in from the top only down to the threshold value region: We recommend that a gravimetric test be performed.

C: The progress bar is filled in to a point well above the threshold value:

The pipette is leaky. A gravimetric test should be performed; if necessary, the unit should be sent in for repair.

What if there are leaks?

If a dynamic test indicates the presence of leaks, you can attempt to localize the source with further testing:

1. A static test should be performed next.
The plunger button is not pressed; the piston is not moved. If the pipette is now leak-tight, this suggests there are scratches on the piston.
2. The dynamic test should be performed without a tip.
Adapter 2 must be used in this case. If the pipette is now leak-tight, it is possible that an unsuitable tip was used.

What do the letters after some volumes mean?

The additional letters after some volumes on the volume selection menu have the following meanings:

- y Yellow (instrument with yellow tip)
- c Clear tip
- n nano-cap tip
- soft During the testing of 5 mL and 10 mL pipettes, the piston might not move all the way back to top in the dynamic test, but might remain stuck in the instrument due to vacuum and friction. If this is the case, select the "5 ml soft or 10 ml soft" test mode.

Note on duration of tests

In the first test after the device has been switched on, a relatively large volume must be pumped out, after which the device must be given time to stabilize (relaxation). Regularly switching between large and small pipette types leads to longer test durations because the measurement pressure must be re-set each time.

The measurement time for pipettes up to 1000 μL is generally 5 s, whereas it can be 20 s for 10 mL instruments.

5.4 Testing multichannel pipettes

Example: 200 µL 8-channel air-displacement pipette, dynamic test

Before performing the test, set the pipette to the nominal volume. Before testing, mount new unused disposable tips and insert four channels with mounted tips into the previously selected and attached Adapter 3. Select an adapter; see Test preparation, p. 58.


What to do	How to do it	Button to use	Display
1. Call up the Main Menu	Press or turn the select knob to call up the Main Menu.		Main menu Pipette type Self-test Device OFF Settings Selection back
2. Set up the self-test	Turn the select knob to the "Pipette Type" setting.		Main menu Pipette type Self-test Device OFF Settings Selection back
3. Confirm the pipette type	Press the "Execute" button. The "Channels" submenu is called up.		Channels Single Channel Multichannel select back
4. Select the Quick Test	Turn the select knob to the "Multichannel" setting.		Channels Single Channel Multichannel select back
5. Confirm the single-channel pipette	Press the "Execute" button. The "Nominal Volume" menu is called up.		Nominal Volume 25 µL 50 µL 100 µL 200 µL select back
6. Select the pipette volume	Turn the select knob to the "200 µL" setting.		Nominal Volume 25 µL 50 µL 100 µL 200 µL select back

What to do	How to do it	Button to use	Display
7. Confirm the pipette volume	Press the "Start/Execute" button. The start submenu is called up.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">ready</p> <p style="text-align: center; margin: 0;">MC- 200µl</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
8. Start the test	Press the "Start/Execute" button. The test is begun, and the test pressure is set (preparing...).		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Test</p> <p style="text-align: center; margin: 0;">preparing...</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
Press the pipette plunger button 2-3 times	Shortly before preparation is finished, the LED light lights up once in red and once in green. Then the progress bar appears and the system automatically starts the measurement. Press the plunger button down evenly 2-3 times during the measurement. The piston must return to the initial position before the measurement is finished! The progress bar indicates the remaining time and the currently achieved vacuum. The LED indicator blinks orange.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Test</p> <div style="background-color: #ccc; width: 100%; height: 15px; margin: 5px 0;"></div> <p style="text-align: center; margin: 0;">-121.21 hPa</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div> <p style="text-align: center; margin-top: 5px;">(pressure value shown as an example)</p> <div style="text-align: center; margin-top: 10px;">  <p style="margin: 0;">LED orange</p> </div>
9. End the test	After the test is finished, the display indicates the leak rate that was determined and whether the system passed the test. To repeat the test, press the "Start/Execute" button. To go back to the Main Menu, turn the select knob.	<div style="text-align: center; margin-bottom: 20px;">  <p style="margin: 0;">LED green</p> </div> <div style="text-align: center;">  <p style="margin: 0;">LED red</p> </div>	<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">PASS</p> <p style="text-align: center; margin: 0;">MC – 200 µl Q_L: 0.04 hPa*ml/s</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> start </div> <p style="text-align: center; margin-top: 5px;">(leak rate shown as an example)</p> <div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">FAIL</p> <p style="text-align: center; margin: 0;">MC – 200 µl Q_L: 4.6 hPa*ml/s</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> start </div> <p style="text-align: center; margin-top: 5px;">(leak rate shown as an example)</p>

NOTICE

You can use the correlation table to estimate the missing volume.

What does the vertical progress bar on the display show?

A mark is added to the progress bar. It indicates the threshold value. It corresponds to approx. $\frac{1}{4}$ of the volume tolerance according to ISO 8655-2.

A: The progress bar is filled in all the way to the bottom: The pipette is leak-tight.

B: The progress bar is filled in from the top only down to the threshold value region: We recommend that a gravimetric test be performed.

C: The progress bar is filled in to a point well above the threshold value:

The pipette is leaky. A gravimetric test should be performed; if necessary, the unit should be sent in for repair.

What if there are leaks?

If a dynamic test indicates the presence of leaks, you can attempt to localize the source with further testing:

1. A static test should be performed next.
The plunger button is not pressed; the piston is not moved.
If the instrument is now leak-tight, this may be an indication of any of the following:
Scratches on the piston
Defective piston seal
Use of unsuitable silicone grease
2. The dynamic test should be performed without a tip.
Adapter 3 can be reused for this.
If the instrument is now leak-tight, possible causes may be any of the following:
Use of unsuitable tips
Defective tip cone
3. If testing a multichannel pipette: Close three channels with the caps provided, and test the individual channels to localize the leak. Select the pipette type "Single Channel" for this!

NOTICE

When using pipettes with pistons lubricated with silicone grease, only use the grease recommended by the manufacturer!

When using pipettes in which the pistons are operated without silicone grease, be sure not to lubricate the pistons; otherwise, irreversible damage may result!
Please refer to the corresponding operating manual.

Note on duration of tests

In the first test after the device has been switched on, a relatively large volume must be pumped out, after which the device must be given time to stabilize (relaxation). Regularly switching between large and small pipette types leads to longer test durations because the measurement pressure must be reset each time.

The measurement time is generally 5 s for pipettes up to 500 μL and 10 s for 1250 μL pipettes.

6 Switching the device off

What to do	How to do it	Button to use	Display
1. Call up the Main Menu	Press the select knob to call up the Main Menu.		<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main menu</p> <p>Pipette type</p> <p>Self-test</p> <p>Device OFF</p> <p>Settings</p> <hr/> <p style="display: flex; justify-content: space-between;">Selection back</p> </div>
2 Select "Switch OFF"	Turn the select knob to the "Switch OFF" setting.		<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main Menu</p> <p>Pipette Type</p> <p>Self-test</p> <p>Switch OFF</p> <p>Settings</p> <hr/> <p style="display: flex; justify-content: space-between;">select back</p> </div>
3. Switch the device off	Press the left "Start/Execute" button to switch the device off.		<p style="font-size: 2em; margin: 0;">OFF</p>
4. Switch the device back on	Press any button to turn the device on again. At first, the display shows the set inspection date for a few seconds. Then it switches automatically to the "Ready" display with the previously tested pipette.		<div style="background-color: #2c3e50; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Inspection date</p> <p style="text-align: center; margin: 0;">1 April 2019</p> <hr/> <p style="text-align: center; margin: 0;">ready</p> <p style="text-align: center; margin: 0;">SC - 100µl</p> <hr/> <p style="display: flex; justify-content: space-between;">select back</p> </div>

7 Make settings

7.1 Select the language

The "Settings" submenu contains the following selectable items: Language, Contrast, Inspection Date, Auto Power Off Time, Pressure Units, and Factory Settings.

What to do	How to do it	Button to use	Display
1. Select "Settings" on the Main Menu	Press the select knob to call up the Main Menu.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main Menu</p> <p>Pipette Type Self Test Switch OFF Settings . . .</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
2. Select "Settings"	Press the "Execute" button. The "Settings" submenu is called up.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
3. Call up the "Language" menu	Turn the select knob to the "Language" setting.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
4. Call up the "Language" menu	Press the "Execute" button. The "Language" menu is called up. Turn the select knob to select a language.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Language</p> <p>Deutsch English Español</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> o.k. back </div>
5. Confirm the language	Press the "Execute" button. The "Settings" menu appears in the language selected.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
	The "Back" button leads to the Main Menu.		


7.2 Set the contrast

The color intensity of the display can be changed interactively.

What to do	How to do it	Button to use	Display
1. Select "Settings" on the Main Menu	Press the select knob to call up the Main Menu.		<div style="background-color: #333; color: #fff; padding: 5px;"> <p style="text-align: center; margin: 0;">Main Menu</p> <p>Pipette Type Self Test Switch OFF Settings . . .</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
2. Select "Settings"	Press the "Execute" button. The "Settings" submenu is called up.		<div style="background-color: #333; color: #fff; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
3. Select the contrast	Turn the select knob to the "Contrast" setting.		<div style="background-color: #333; color: #fff; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
4. Call up the contrast	Press the "Execute" button. The "Contrast" menu is called up.		<div style="background-color: #333; color: #fff; padding: 5px;"> <p style="text-align: center; margin: 0;">Contrast</p> <div style="display: flex; align-items: center; justify-content: center;">  50 </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> o.k. back </div>
5. Set the contrast	Set the contrast to a value from "10" to "100" by turning the select knob.		<div style="background-color: #333; color: #fff; padding: 5px;"> <p style="text-align: center; margin: 0;">Contrast</p> <div style="display: flex; align-items: center; justify-content: center;">  76 </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> o.k. back </div>
6. Confirm the contrast	Press the "Execute" button. The "Settings" menu reappears.		<div style="background-color: #333; color: #fff; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> select back </div>
	The "Back" button leads to the Main Menu.		

7.3 Inspection date

The date on when the PLT unit should be sent to BRAND or an authorized repair shop for testing can be set individually. This date is displayed each time the unit is switched on.

What to do	How to do it	Button to use	Display
1. Select "Settings" on the Main Menu	Press the select knob to call up the Main Menu.		Main Menu Pipette Type Self Test Switch OFF Settings . . . select back
2. Confirm the settings	Press the "Execute" button. The "Settings" submenu is called up.		Settings Language Contrast Inspection date Auto-Power-Off Time select back
3. Select the inspection date	Turn the select knob to the "Inspection Date" setting.		Settings Language Contrast Inspection date Auto-Power-Off Time select back
4. Call up the inspection date	Press the "Execute" button. The previously set inspection date is displayed.		Inspection date 23 January 2023 o.k. back
5. Set the inspection date	Press the select knob to switch between year, month, and day. The selected position is shown with a frame around it. (To go back to the "Settings" menu, press the select knob with the "Day" position selected.) Turn the select knob to change the value of the selected parameter.		Inspection date 23. January 2023 o.k. back
6. Confirm the inspection date	Press the "Execute" button. The display returns automatically to the "Settings" menu.		Settings Language Contrast Inspection date Auto-Power-Off Time select back

What to do	How to do it	Button to use	Display
	To go back to the Main Menu, press the "Back" button.		

7.4 Auto Power Off Time

The auto power off time of the PLT unit can be set individually in minutes ranging from "OFF" to 30 min.

What to do	How to do it	Button to use	Display												
1. Select "Settings" on the Main Menu	Press the select knob to call up the Main Menu.		<table border="1"> <thead> <tr> <th colspan="2">Main Menu</th> </tr> </thead> <tbody> <tr> <td colspan="2">Pipette Type</td> </tr> <tr> <td colspan="2">Self Test</td> </tr> <tr> <td colspan="2">Switch OFF</td> </tr> <tr> <td colspan="2">Settings . . .</td> </tr> <tr> <td>select</td> <td>back</td> </tr> </tbody> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings . . .		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings . . .															
select	back														
2. Confirm the settings	Press the "Execute" button. The "Settings" submenu is called up.		<table border="1"> <thead> <tr> <th colspan="2">Settings</th> </tr> </thead> <tbody> <tr> <td colspan="2">Language</td> </tr> <tr> <td colspan="2">Contrast</td> </tr> <tr> <td colspan="2">Inspection date</td> </tr> <tr> <td colspan="2">Auto-Power-Off Time</td> </tr> <tr> <td>select</td> <td>back</td> </tr> </tbody> </table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		select	back
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
select	back														
3. Select the auto power off time	Turn the select knob to the "Auto Power Off Time" setting.		<table border="1"> <thead> <tr> <th colspan="2">Settings</th> </tr> </thead> <tbody> <tr> <td colspan="2">Language</td> </tr> <tr> <td colspan="2">Contrast</td> </tr> <tr> <td colspan="2">Inspection date</td> </tr> <tr> <td colspan="2">Auto-Power-Off Time</td> </tr> <tr> <td>select</td> <td>back</td> </tr> </tbody> </table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		select	back
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
select	back														
4. Call up the auto power off time	Press the "Execute" button. The previously set auto power off time is displayed.		<table border="1"> <thead> <tr> <th colspan="2">Auto-Power-Off Time</th> </tr> </thead> <tbody> <tr> <td colspan="2">3 min</td> </tr> <tr> <td>o.k.</td> <td>back</td> </tr> </tbody> </table>	Auto-Power-Off Time		3 min		o.k.	back						
Auto-Power-Off Time															
3 min															
o.k.	back														
5. Set the inspection date	Turn the select knob. The auto power off time can be set to a value ranging from "OFF" to 30 min.		<table border="1"> <thead> <tr> <th colspan="2">Auto-Power-Off Time</th> </tr> </thead> <tbody> <tr> <td colspan="2">15 min</td> </tr> <tr> <td>o.k.</td> <td>back</td> </tr> </tbody> </table>	Auto-Power-Off Time		15 min		o.k.	back						
Auto-Power-Off Time															
15 min															
o.k.	back														

What to do	How to do it	Button to use	Display
6. Confirm the inspection date	Press the "Execute" button. The display returns automatically to the "Settings" menu.		Settings Language Contrast Inspection date Auto-Power-Off Time select back
	The "Back" button leads to the Main Menu.		

7.5 Pressure units (hPa, mbar, torr)

What to do	How to do it	Button to use	Display
1. Select "Settings" on the Main Menu	Press the select knob to call up the Main Menu.		Main Menu Pipette Type Self Test Switch OFF Settings . . . select back
2. Confirm the settings	Press the "Execute" button. The "Settings" submenu is called up.		Settings Contrast Inspection date Auto-Power-Off Time Pressure unit select back
3. Select the pressure units	Turn the select knob to the "Pressure Units" setting.		Settings Inspection date Auto-Power-Off Time Pressure unit Factory settings select back
4. Call up the pressure units	Press the left "Execute" button. The currently set pressure units appear.		Pressure unit mbar o.k. back
5. Set the pressure units	Turn the select knob to select hPa, mbar, or torr as the pressure units.		Pressure unit hPa o.k. back

What to do	How to do it	Button to use	Display
6. Confirm the pressure units	Press the "Execute" button. The display returns automatically to the "Settings" menu.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure unit Factory settings</p> <div style="border: 1px solid white; width: 100%; height: 15px; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> select back </div> </div>
	The "Back" button leads to the Main Menu.		

7.6 Factory settings

What to do	How to do it	Button to use	Display
1. Select "Settings" on the Main Menu	Press the select knob to call up the Main Menu.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Main Menu</p> <p>Pipette Type Self Test Switch OFF Settings . . .</p> <div style="border: 1px solid white; width: 100%; height: 15px; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> select back </div> </div>
2. Confirm the settings	Press the "Execute" button. The "Settings" submenu is called up.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure unit Factory settings</p> <div style="border: 1px solid white; width: 100%; height: 15px; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> select back </div> </div>
3. Select the factory settings	Turn the select knob to the "Factory Settings" setting.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure unit Factory settings</p> <div style="border: 1px solid white; width: 100%; height: 15px; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> select </div> </div>
4. Call up the factory settings	Press the "Execute" button. The confirmation prompt appears.		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Factory settings</p> <p style="text-align: center; margin: 10px 0 10px 40px;">Please confirm</p> <div style="display: flex; justify-content: space-between; font-size: small;"> o.k. back </div> </div>
5. Confirm the factory settings	Press the "Execute" button. The display changes to the "English" language selection. After the language has been selected and confirmed, the		<div style="background-color: #333; color: white; padding: 5px;"> <p style="text-align: center; margin: 0;">Language</p> <p>Deutsch English Español</p> <div style="border: 1px solid white; width: 100%; height: 15px; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> o.k. back </div> </div>

What to do	How to do it	Button to use	Display
	display automatically changes back to the "Settings" menu.		<div data-bbox="804 113 1025 248" style="background-color: #333; color: white; padding: 5px;"> <p>Settings</p> <p>Inspection date</p> <p>Auto-Power-Off Time</p> <p>Pressure unit</p> <p style="border: 1px solid white; display: inline-block; padding: 2px;">Factory settings</p> </div> <div data-bbox="804 248 1025 288" style="background-color: #333; color: white; padding: 2px;"> <p>select</p> </div>
	The "Back" button leads to the Main Menu.		

8 Firmware update

Firmware updates, e.g., due to inclusion of additional pipette types (new volumes), are made available by BRAND on CD.

Registered users are informed automatically.

Unregistered users receive this update on request.

Installation:

- a. Place the CD-ROM in the drive (drive "D" is used in this example).
- b. Click on "Start" on the desktop, and enter "D:\setup.exe" in the "Execute" command field to install the firmware update.
- c. Follow the instructions on the screen.

USB port – see Operating elements, p. 52

9 Correlation table

Relationship between QL and volume loss (approx.)

For a leak of a given size, the volume loss in the pipette during pipetting depends on various factors, e.g., the tip opening, and can thus vary significantly. (Tolerance for the measured leak rate: $\pm 20\%$)

Correlation table

Pipette tested			Leak rate QL		Pipette volume loss (approx.)
Channels	Nominal Volume		[hPa*m/s]	[μ l]	
1-channel	1	μ L	4		0.05
1-channel	2	μ L	4		0.09
1-channel	2.5	μ L	6		0.13
1-channel	5	μ l y	12		0.1
1-channel	5	μ L	4		0.1
1-channel	10	μ l y	4		0.1
1-channel	10	μ l c/n	3		0.1
1-channel	20	μ l y	19		0.2
1-channel	20	μ l c	6		0.2
1-channel	25	μ L	16		0.5
1-channel	50	μ L	10		0.5

Pipette tested			Leak rate QL		Pipette volume loss (approx.)
1-channel	100	μ L	14		0.8
1-channel	200	μ L	13		2
1-channel	250	μ L	13		4
1-channel	300	μ L	12		4
1-channel	500	μ L	10		4
1-channel	1000	μ L	12		8
1-channel	1250	μ L	26		16
1-channel	1500	μ L	36		16
1-channel	2000	μ L	40		16
1-channel	2500	μ L	66		40
1-channel	5	ml	52		40
1-channel	5	ml soft	21		41
1-channel	10	ml	39		60
1-channel	10	ml soft	20		61
4-channel	1	μ L	18		0.1
4-channel	2	μ L	19		0.2
4-channel	2.5	μ L	29		0.2
4-channel	5	μ L	22		0.2
4-channel	10	μ L	8		0.2
4-channel	20	μ L	10		0.4
4-channel	25	μ L	21		1
4-channel	50	μ L	13		1
4-channel	100	μ L	13		2
4-channel	200	μ L	14		3
4-channel	250	μ L	44		8
4-channel	300	μ L	47		8
4-channel	500	μ L	29		8
4-channel	1000	μ L	25		16
4-channel	1250	μ L	25		16

10 Troubleshooting


Errors occurring during pipette tests

Fault	Error display	Possible cause	Remedy	
Test interrupted	<div style="background-color: #333; color: white; padding: 5px; text-align: center;"> FAIL Test aborted Leak rate cannot be measured! </div>	No pipette in adapter or pipette not inserted firmly enough	Insert pipette more firmly.	
		Tested pipette full of leaks	Pipette is defective and cannot be tested with this testing unit. Send pipette in for repair.	
		back		
		Tip not completely sealing shaft	Mount another tip. Check shaft for damage.	
		Wrong pipette type set	Select correct pipette type from menu.	
		Pipette adapter not leak-tight	Check adapter for lint or abrasion and check filter for moisture. Clean adapter and replace filter. Perform self-test.	
		Pipette adapter loose	Tighten lock nut. Perform self-test	
		Moisture in measurement room; adapter filter clogged	Replace adapter filter; perform self-test (Quick Test).	
		For 5 mL and 10 mL pipettes: Measurement pressure too high for pipette's piston spring Piston being pulled downwards during evacuation by pump	Set pipette type to "5 ml soft" or "10 ml soft".	
Measured leak rate implausible		Moisture or leak in PLT unit	Perform self-test (Full Test).	
PLT not responding at all or behaving unusually		Software error; electrical fault	Disconnect all cables from device. Reconnect power cable after approx. 5 seconds.	

NOTICE

If errors cannot be resolved through the specified measures, send the device in for repair.

Errors occurring during self-test

Error code	Cause	Remedy
1000 1010	Pressure sensor defective	Send device in for repair.
1020	Pump running continuously	Send device in for repair.
1040 1050	Pumping out taking too long	Check whether cap is securely inserted and adapter is screwed in tightly. If so, then pump is defective. Send device in for repair.
1030 1060 1070 1080 1090	Valve defective	Send device in for repair.
1020	Aeration through aeration channel taking too long	Send device in for repair.
1150 1160	Device leaking	Check cap, adapter seating, and adapter filter, and replace if necessary. Repeat self-test. If errors cannot be resolved, send the device in for repair.
1170	Aeration through adapter taking too long	Remove cap if not already done. If adapter filter is contaminated, replace filter.
1200	Electronics fault	Send device in for repair.

11 Product markings

Symbol or Number	Meaning
	With this mark, we confirm that the product complies with the requirements set out in the EC Directives and has undergone testing according to the specified test procedures.
	UKCA: United Kingdom Conformity Assessed With this mark, we confirm that the product complies with the requirements set out in the UK Designated Standards.
www.brand.de/ip	Patent information
XXNXXXXX	Serial number
	Heed the information on the device, on the accessories, and in the operating manual.
	China RoHS (EFUP) EFUP defines the period in years during which the hazardous substances contained in electrical and electronic equipment will not leak or mutate under normal operating conditions. When used normally by the user, such electrical and electronic products will not lead to serious environmental pollution, serious bodily injury, or damage to the user's physical assets.
	The electrical device should not be disposed of with household waste.

12 Ordering information

PLT unit (Pipet Leak Testing Unit)

Description	Order No.
PLT unit (pipette leak testing unit), including 1 adapter each for testing air-displacement pipettes with tips (mounted) and without tips, 2 caps, 3 replacement PE filters for the pipette adapters, 1 universal AC adapter, and 1 quality certificate	703970

Pipette adapter for testing with tips

Description	Order No.
Pipette adapter for testing single-channel pipettes with and without tips, including 1 cap	703975

Pipette adapter for testing without tips

Description	Order No.
Pipette adapter for testing single-channel pipettes without tips, including 1 cap	703976

4-channel pipette adapter

Description	Order No.
4-channel pipette adapter for testing multichannel pipettes with and without tips, including 4 caps	703977

Filters for the pipette adapters

Description	Order No.
PE filters + cap for pipette adapters (10 filters + 1 cap)	703978

Universal AC adapter

Description	Order No.
Universal AC adapter Input: AC 100 V - 240 V, 50/60 Hz Output: DC 6.5 V, 800 mA	703979

PLTconnect

Description	Order No.
PLTconnect software For the documentation of tests (pack of 1)	703980

13 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The complete SOP for testing can be downloaded from www.brand.de or www.brandtech.com.

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at www.brand.de (Service & Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

14 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

15 Repairs

15.1 Sending for repair

NOTICE

Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Outside USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791


China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn

info@brand.co.in
www.brand.co.in

www.brand.cn.com

16 Disposal

This symbol means that at the end of their service life, batteries/accumulators and electronic devices must be disposed of separately from household waste (unsorted municipal waste).

Electronic devices must be disposed of in accordance with Directive 2012/19/ EU of the European Parliament and of the Council from July 04, 2012 on waste from electrical and electronic equipment and in compliance with national disposal regulations.

Both batteries and accumulators (rechargeable batteries) contain materials that can be damaging to the environment and human health. Therefore, they must be properly disposed of in accordance with Directive 2006/66/EC of the European Parliament and of the Council from September 06, 2006 on batteries and accumulators and in compliance with national disposal regulations. Only dispose of fully discharged batteries and accumulators.

Subject to technical changes, errors, and misprints.

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