Gebrauchsanleitung | Operating manual | Mode
d'emploi | Instrucciones de manejo | Istruzione | In-
struções | 操作手册

Dispensette® S Trace Analysis

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Indien:

BRAND Scientific Equipment Pvt. Ltd.
 303, 3rd Floor, 'C' Wing, Delphi
 Hiranandani Business Park,
 Powai
 Mumbai-400 076 (India)
 T +91 22 42957790
 F +91 22 42957791
 info@brand.co.in
 www.brand.co.in

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Table of contents

1 Introduction	33	14.1 Sending for repair.....	58
1.1 Scope of delivery.....	33	15 Warranty	60
1.2 Terms of use.....	33	16 Disposal	60
2 Safety Instructions	34		
2.1 General safety instructions.....	34		
2.2 Function.....	35		
2.3 Limitations of Use	35		
2.4 Operating Limitations.....	35		
2.5 Operating Exclusions	35		
2.6 Storage Conditions	36		
2.7 Recommended application range...	36		
3 Functional and operating elements	37		
4 Assembly	38		
4.1 First Steps.....	38		
4.2 Priming	40		
5 Operation	42		
5.1 Dispensing.....	42		
5.2 Accessories	43		
5.3 Bottle stand	46		
6 Error limits	46		
7 Checking the Volume (Calibration)	47		
8 Adjustment	48		
8.1 Making adjustments.....	48		
8.2 Adjustment range.....	48		
9 Cleaning	49		
9.1 Standard cleaning.....	49		
9.2 Cleaning for trace analysis.....	50		
9.3 Replacing discharge tube/valves	50		
9.4 Replacing the dispensing cartridge.	52		
10 Troubleshooting	53		
11 Product markings	54		
12 Ordering Information	55		
13 Accessories/spare parts	55		
14 Repairs	58		

1 Introduction

1.1 Scope of delivery

Bottle-top dispenser Dispensette® S Trace Analysis, telescopic filling tube, discharge tube or discharge tube with recirculation valve and recirculation tube (included only in recirculation valve models), assembly tool, three bottle adapters, a performance certificate and this operating manual.

Nominal volume [ml]	Adapter for bottle thread	Filling tube Length [mm]
10	GL 28/ S 28 (ETFE), GL 32 (ETFE), GL 38 (ETFE), S 40 (PTFE)	125-240

1.2 Terms of use

- Please carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: www.brand.de.

1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2.2 Symbols

Symbol	Meaning
	Hazardous area

1.2.3 Format

Layout	Meaning	Layout	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇨	Indicates a result.

2 Safety Instructions

2.1 General safety instructions

Please read carefully!

The can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual before operation.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
3. Observe all specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
5. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations. Comply with the operating exclusions; see Operating Exclusions, p. 35. If in doubt, contact the manufacturer or supplier.
6. Always perform work in a manner that does not endanger yourself or other people. When dispensing, the discharge tube must always point away from you or any other person. Avoid splattering. Only use suitable vessels.
7. Never press down the piston when the discharge tube closure is attached.
8. Never remove the discharge tube while the dispensing cylinder is filled.
9. Reagents can accumulate in the screw cap of the discharge tube. Thus, the screw cap should be cleaned regularly.
10. To prevent tipping, use a bottle stand – particularly with small bottles and when using the flexible discharge tube.
11. When mounted to a bottle, never carry the instrument by the cylinder sleeve or the valve block. Breakage or loosening of the cylinder can lead to personal injury from chemicals, see First Steps, p. 38, Fig. 4.
12. Never use force. Use smooth gentle movements to operate the piston upwards and downwards.
13. Use only original accessories and original replacement parts. Do not make any technical modifications. Do not dismantle the instrument any further than is described in the operating manual!
14. Always check that the instrument is in proper working condition before use. The user can come into contact with media if the instrument has been insufficiently cleaned or inspected. If there is a sign of a potential malfunction (e.g., piston difficult to move, sticking valves or leakage), immediately stop dispensing and consult the Troubleshooting, p. 53 section. Contact the manufacturer, if necessary. In case of discolorations, check if any material fatigue is present. If there is any doubt, replace components.
15. Always tighten the retaining ring between the valve block and dispensing cylinder by hand. Do not use any tools.
16. Autoclaving the instrument is not permitted!

2.2 Function

With the Dispensette® S Trace Analysis bottle-top dispenser, liquids can be dispensed directly from the supply bottle. The instruments are marked DE-M and optionally equipped with recirculation valve.

2.2.1 Operation

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

Various fluoroplastics (e.g. ETFE, FEP, PFA, PCTFE, PTFE), Al₂O₃-sapphire, platinum-iridium or tantalum depending on the version (see label on the dispensing tube and the recirculation valve).

2.3 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- Operating temperature from +15 °C to +40 °C (from 59 °F to 104 °F) of instrument and reagent
- Vapor pressure up to max. 600 mbar. Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling
- Kinematic viscosity up to 500 mm²/s (dynamic viscosity [mPas] = kinematic viscosity [mm²/s] x density [g/cm³])
- Density up to 3.8 g/cm³

2.4 Operating Limitations

- Liquids, which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or concentrated alkaline solutions). If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately (Cleaning, p. 49).
- When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
- The instrument is designed for general laboratory applications and complies with the relevant standards, e.g. DIN EN ISO 8655. Compatibility of the instrument for a specific application (e.g., trace material analysis, food sector etc.) must be checked by the user. Approvals for specific applications, e.g. for production and administration of food, pharmaceuticals or cosmetics are not available.

2.5 Operating Exclusions

2.5.1 Dispensette® S Trace Analysis

Never use Dispensette® S Trace Analysis for:

- Liquids attacking Al₂O₃-sapphire or fluoroplastics such as ETFE, FEP, PFA PCTEF and PTFE (e.g., dissolved sodium azide*)
- organic solvents
- Trifluoroacetic acid
- Explosive liquids (e.g., carbon disulfide)

- Liquids which are decomposed catalytically by platinum-iridium (e.g., H_2O_2) or attack tantalum. Take notice of the valve spring material of the instrument version in use.
- Suspensions (e.g., of charcoal), as solid particles may clog or damage the instrument

* Dissolved sodium azide is permitted up to a concentration of max. 0.1 %.

2.6 Storage Conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

2.7 Recommended application range

Select valves with valve springs made of platinum-iridium or tantalum according to the intended use. The instrument can be used for the following dispensing media:

Dispensing medium	Valve springs: Pt-Ir	Valve springs: Ta
Ammonia solution	+	+
Bromine	+	+
Acetic acid	+	+
Hydrofluoric acid*	+	-
Sodium hydroxide, 30 %	+	-
Perchloric acid	+	+
Phosphoric acid	+	+
Nitric acid	+	+
Hydrochloric acid	+	+
Sulphuric acid	+	+
Water	+	+
Hydrogen peroxide	-	+

+ suitable - not suitable

* Hydrofluoric acid slightly attacks the sapphire components (99.99% Al_2O_3) used in the valves and can release aluminum ions. Depending on the analysis method used (inorganic trace analysis), this can lead to slightly increased aluminum blank values. To reduce the aluminum values, we recommend discarding 3 to 5 doses of 2 ml each before the analysis.

Compounds that contain fluoride, such as NaF, attack tantalum.

This table has been carefully tested and is based on the most current information available. Always observe the operating manual of the instrument and the specifications provided by the reagent manufacturer. If you need chemical hazard statements that are not on the list, you are welcome to contact BRAND.

Last updated: 0419/2

3 Functional and operating elements

- 1 Cover plate
- 2 Piston seat
- 3 Housing
- 4 Pointer
- 5 Volume adjustment wheel
- 6 Retaining ring
- 7 Recirculation valve (optional)
- 8 Valve block
- 9 Valve block adapter (GL 45 bottle thread)
- 10 Telescoping filling tube
- 11 Recirculation tube (optional)
- 12 Screw cap
- 13 Discharge tube
- 14 Valve lever, recirculation valve
- 15 Dosing unit with serial number

 The retaining ring between the valve block and the dosing unit must always be securely tightened.

Telescoping filling tube and recirculation tube

Assembly tool

4 Assembly

4.1 First Steps

⚠ WARNING

Follow the safety instructions

- > Wear protective clothing, eye protection and protective gloves!
- > Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- > Follow all safety instructions and comply with the limitations of use, see Limitations of Use, p. 35.
- > Comply with the operating limitations, see Operating limitations, p. 35.

NOTICE

Selecting the correct discharge valve and discharge tube

The discharge valve and discharge tube are labeled with the spring material. The label 'Pt-Ir' or 'Ta' must be the same on each instrument. The spring material derives from the application range (observe Recommended application range, p. 36).

When labels do not match, the desired application is not possible as components can become corroded or decomposed.


1. Checking the retaining ring

a.

- a. Check that the retaining ring is screwed tight.

2. Mounting the filling tube/recirculation tube

1	Opening for recirculation tube
2	Olive for filling tube

- a. Adjust the length of the telescoping filling tube to the bottle height and attach it.
- b. Center and attach the filling tube (part with smaller diameter) carefully to avoid damaging the nozzle.
 - ⇒ If a discharge tube with a recirculation valve is used, the recirculation tube must also be installed.
- c. Insert the recirculation tube with the opening pointing outward.

3. Mounting and aligning the instrument on a bottle

NOTICE

Cleaning the instrument before use in trace analysis

Before use in trace analysis, the instrument must be thoroughly cleaned. See *Cleaning for trace analysis*, p. 50.

NOTICE

Prevent tipping

To prevent tipping, use a bottle stand – particularly with small bottles and when using the flexible discharge tube.

- a. Screw the instrument (GL 45 threads) onto the reagent bottle, and then align the discharge tube with the bottle label. This is done by rotating the valve block with the discharge tube.

4. Transporting the instrument

⚠ WARNING

Wear protective gloves

Always wear protective gloves when touching the instrument or the bottle, especially when using hazardous media (e.g., HF).

Incorrect handling

Incorrect handling can, among other things, lead to breakage of the bottleneck, retaining ring or the bottle adapter.

- a. For bottles with other thread sizes, select a suitable bottle adapter.
- b. Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- c. When mounted to a reagent bottle, always carry the instrument as shown in the figure!

4.2 Priming

⚠ WARNING

Follow during every use, especially with hazardous media

- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 35 and Operating exclusions, p. 50.

NOTICE

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. To avoid splashes dispense slowly. Carry out cleaning for trace analysis, depending on requirements (Cleaning for trace analysis, p. 50).

Instruments with recirculation valve

- a. Open the screw cap of the dispensing tube. To avoid splashes, hold the discharge tube orifice on the inner wall of a suitable receiving vessel.
- b. Turn valve to 'Recirculate'.
- c. For deaerating gently pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this process at least 5 times.
- d. Turn valve to 'Dispense'.
- e. To avoid splashes when deaerating hold the discharge tube on the inner wall of a suitable receiving vessel and dispense liquid to vent deaerate the discharge tube until it is bubble-free. Wipe away any remaining drops from the discharge tube.


Instruments without recirculation valve

- a. Open the screw cap of the discharge tube (see 'instrument with recirculation valve', Fig. a.). To avoid splashes, hold discharge tube orifice on the inner wall of a suitable receiving vessel.
- b. For deaerating pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure approximately 5 times until the discharge tube is bubble-free.

5 Operation

5.1 Dispensing

1. Selecting the volume

- a. Loosen the volume selector thumb screw $\frac{3}{4}$ turn (1), set the pointer to the desired volume (2) and then retighten the volume thumb screw (3).

2. Dispensing

⚠ WARNING

Follow during every use, especially with hazardous media


- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 35 and Operating exclusions, p. 35.

- a. Unscrew the screw cap of the dispensing tube.
- b. When using instruments equipped with a recirculation valve, turn the valve to 'Dispensing'.
- c. Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.

- d. Gently lift the piston until the upper stop and then depress piston slowly and steadily with minimal force until the lower stop.

- e. Wipe off the discharge tube against the inner wall of the receiving vessel.
- f. Seal the dispensing tube with the screw cap.

NOTICE

After using the piston, always press it down to the lower stop (parking position). If the piston is not pressed down to the lower stop, unintentional media leaks may occur.

NOTICE

The filled status of the instrument must be specially marked during cleaning.

5.2 Accessories

5.2.1 Flexible discharge tube with recirculation valve

With the exception of HF, the flexible discharge tube can be used for serial dispensing (Accessories/ spare parts, p. 55).

The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. Before use, the tube must lie in regular loops and must not be twisted. The operating exclusions for the respective instrument in use apply.

Assembly

⚠ WARNING

Use only undamaged tubing

There should be no visible damage to the discharge tube (e.g. kinks, etc.). Each time you are going to use the tubing, examine it carefully.




- > To dispense aggressive liquids, you should take safety measures in addition to the normal precautions. We recommend use of a protective shield.
- > The bottle must be supported using a bottle stand.
- > To help avoid reagent splashing from the tube, always grip the tube firmly by the handle and replace into the holder after use.
- > For cleaning rinse the tube carefully.
- > Do not dismantle!

Unsuitable for hydrofluoric acid (HF)

The flexible discharge tube must never be used to dispense HF (hydrofluoric acid)!

Checking for the correct valve marking

Always take note of the valve marking! (See Replacing the discharge valve, p. 51).

Prerequisite:

- If the instrument was in use, it must be cleaned before installing the flexible discharge tube (Cleaning, p. 49).
- a. Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.
- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.

- c. Push the flexible discharge tube holder from the bottom of the valve block and tighten it. In doing so, the instrument must not be mounted to a bottle. Install the receiver tube.
- d. Press the plug of the recirculation valve downwards.

- e. Slide the flexible discharge tube housing into the valve block up to the stop.

- f. Push the housing all the way down.
- g. Attach the valve handle that fits the discharge valve and press it in firmly. Note the color coding and marking.

NOTICE

Use a bottle stand (Accessories/spare parts, p. 55).

5.2.2 Drying tube

Use of a drying tube, filled with a suitable absorbent (purchased separately), might be necessary for moisture- and CO₂- sensitive media.

(Accessories/spare parts, p. 55)

Assembly

- a. Use a coin to unscrew the air vent cap.

- b. Screw in the filled drying tube.

- c. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

NOTICE

If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

5.2.3 Sealing ring for valve block

For highly volatile media we recommend to seal the connection from valve block to bottle with the PTFE sealing ring and PTFE tape (Accessories/spare parts, p. 55).

Assembly

- a. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

5.3 Bottle stand

For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over (Accessories/spare parts, p. 55).

Assembly

- a. Position the mounting plate at the appropriate height.
- b. As shown, insert the instrument firmly into the holder until it audibly clicks into place.
- c. Then lock the holder into place with the screw.

6 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. Testing takes place according DIN EN ISO 8655-6 with a completely filled instrument and with uniform and smooth dispensing.

Error limits

Nominal volume ml	A* ≤ ± %	μl	CV* ≤ %	μl
10	0.5	50	0.1	10

*A = Accuracy, CV = Coefficient of Variation

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

Partial volumes

The percentage values for A and CV are relative to the nominal volume (V_N) and must be converted for partial volumes (V_p).

e.g.	Volume	A* ≤ ± %	μl	CV* ≤ %	μl
V_N	10.0	0.5	50	0.1	10
$V_T = 50\% N$	5.0	1.0	50	0.2	10
$V_T = 10\% N$	1.0	5.0	50	1.0	10

*A = Accuracy, CV = Coefficient of Variation

NOTICE

The error limits in DIN EN ISO 8655- 5 are significantly lower. The maximum error for a single measurement can be approximated from the sum of error limits $EL = A + 2 \times CV$ (for 10 ml: $50 \mu\text{l} + 2 \times 10 \mu\text{l} = 70 \mu\text{l}$).

7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to correspond with individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de.

For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from www.brand.de.

Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Error limits, p. 46) is performed as follows:

1. Prepare the instrument

Clean the instrument (Cleaning, p. 49), fill it with distilled H₂O and then prime it carefully.

2. Check the volume

- 10 dispensing operations with distilled H₂O in 3 volume ranges (100 %, 50 %, 10 %) are recommended
- For discharge, depress the piston slowly and steadily without force until the lower stop
- Wipe off the tip of discharge tube.
- Weigh the dispensed quantity on an analytical balance. (Please refer to the user manual of the scale manufacturer.)
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

3. Calculation


x_i = Weighing results	n = Number of weighings
Z = Correction factor (e.g. 1.0029 µl/mg at 20 °C, 1013 hPa)	
Mean	$\bar{x} = \frac{\sum x_i}{n}$
Mean volume	$\bar{V} = \bar{x} \cdot z$
Accuracy*	$A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$
	V_0 = Nominal volume
Coefficient of variation* CV%	$= \frac{100 s}{\bar{V}}$
Standard deviation*	$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$

8 Adjustment


After a long period of usage, an adjustment of the instrument might be necessary.

- Calibrate for example at nominal volume (Checking the Volume (Calibration), p. 47).
- Calculate mean volume (actual value) (Checking the Volume (Calibration), p. 47).
- Adjust the instrument (to the calculated mean volume).
- After the adjustment, further calibration is necessary to confirm the appropriate adjustment.

8.1 Making adjustments

- Insert the pin of the assembly tool into the cover plate (Pos. 1), and break it off with a rotating motion. Discard the adjustment cover.

- Insert the pin of the assembly tool into the adjustment screw and rotate to the left in order to increase the dispensing volume, or rotate to the right to decrease the dispensing volume (e.g. for an actual value of 9.97 ml, rotate approx. 1/2 turn to the left).

- Adjustment is completed.
- ⇒ The change in the adjustment is indicated by a red disk (circle in figure).

8.2 Adjustment range

Max $\pm 60 \mu\text{l}$

One rotation corresponds to $\sim 80 \mu\text{l}$.

9 Cleaning

⚠ WARNING

Components filled with reagent

The cylinder, valves, telescopic filling tube and discharge tube contain reagent!

- > Never remove the discharge tube while the dispensing cylinder is filled.
- > Point the valves and tube openings away from your body.
- > Wear protective clothing, eye protection and protective gloves!

The instrument must be cleaned in the following situations to assure correct operation:

- prior to first use
- immediately when the piston is difficult to move
- prior to a reagent change
- prior to long term storage
- prior to dismantling the instrument
- prior to valve replacement
- regularly when using liquids which form deposits (e.g., crystallizing liquids)
- regularly when liquids accumulate in the screw cap

Autoclaving the instrument is **not** permitted!

9.1 Standard cleaning

1. Completely emptying the instrument

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

2. Rinsing the instrument

- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

3. Rinsing in recirculating position (optional)

If the instrument is equipped with a recirculation valve, it must also be rinsed in recirculating position after rinsing the instrument.

- a. Set the recirculation valve to "recirculate" and fill and discharge the instrument several times.

9.2 Cleaning for trace analysis

Before use in trace analysis, the instrument must first be thoroughly cleaned. In doing so, use reagents having the purity level "per analysis" or better. If contamination of the bottle contents should be avoided, do not use the instrument with a recirculation valve. If the instrument is used with a recirculation valve, the cleaning must be done using both the dispensing and recirculating function (see Figures 1 and 2).

Good results have been achieved in practice with the following recommended cleaning methods. Modify accordingly, as necessary.

Dosing function

Recirculation function

- a. Screw the device onto a bottle filled with **acetone**, deaerate and fill up to the maximum. Leave the piston at the upper stop and close the discharge tube with the screw cap. Dispense twice after approx. 24 hours reaction time, then the device must be completely drained and rinsed 5 times with pure water.
- b. Screw the device onto a bottle filled with approx. **20 % hydrochloric acid**, deaerate and fill up to the maximum. Leave the piston at the upper stop and close the discharge tube with the screw cap.
- c. After approx. 24 hours reaction time, dispense twice and fill again to the maximum.
- d. Repeat step 3 two more times. After another 24 hours of reaction time, completely drain the device and rinse it 5 times with pure water.
- e. Repeat steps 2 to 4 with a solution of approx. **30% nitric acid**.
- f. Screw the device onto a bottle filled with the desired dispensing medium, deaerate and fill up to the maximum. Leave the piston at the upper stop and close the discharge tube with the screw cap.
- g. After approx. 24 hours reaction time, dispense twice and fill again to the maximum.
- h. Repeat step 7 two more times. After an additional 24 hours reaction time, dispense twice and leave the piston at the lower stop.

NOTICE

If the cleaning is insufficient, repeat the cleaning process.

9.3 Replacing discharge tube/valves

⚠ WARNING

The valves, telescoping filling tube and discharge tube contain reagent!

Contact with potentially hazardous media is possible.

- > Clean the instrument before the discharge tube, valves or dosing unit are exchanged.
- > Do not disassemble the dosing unit.
- > Follow the safety instructions (Safety Instructions, p. 34).


NOTICE

A functional test must always be carried out after exchanging components.

9.3.1 Replacing discharge tube

- a.** Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.

- b.** Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.
- c.** Hold coupling piece of the new discharge tube and pull housing up. Push housing into the valve block until it meets the stop.

- d.** Slide the discharge tube housing all the way down.
- e.** For instruments with a recirculation valve, pull up the valve lever to the 'Recirculate' position, and press it in tightly.

9.3.2 Replacing valves


9.3.2.1 Replacing the discharge valve

NOTICE

Always install the valve intended for the particular instrument model and size.

The filling valve of the is the same on all instrument versions; however, the discharge valve is different. Make sure to use only the filling valve marked with 'S' (sapphire). For easy identification, the discharge valves for the are marked with 'Pt-Ir' or 'Ta'.

a.

- a. After disassembling the discharge tube (see Replacing discharge tube, p. 51), use the assembly tool to unscrew the discharge valve.
- b. Firmly screw in the new discharge valve first by hand and then tighten it with the assembly tool. The thread should not be visible.

9.3.2.2 Replacing the filling valve

b.

- a. Pull out the recirculation tube and the telescopic filling tube.
- b. Use the assembly tool to unscrew the filling valve.
- c. Screw in the new filling valve first by hand and then tighten it with the assembly tool.

9.3.2.3 Dislodging a stick valve ball

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200 μ l plastic pipette tip (see the figure at the side).

9.4 Replacing the dispensing cartridge

The dispensing cartridge is a wear part and must be replaced according to the specific conditions of use. The replacement interval depends, for example, on the medium used and number of dispensing cycles, whereby replacement is to be expected after 10,000 full strokes (medium: distilled water). If drops of media are observed on the outside of the dispensing cylinder or in the retaining ring, immediately check the instrument and, in case of doubt, replace the dispensing cartridge.

⚠ WARNING**The valves, telescoping filling tube and discharge tube contain reagent!**

Contact with potentially hazardous media is possible.

- Clean the instrument before the discharge tube, valves or dosing unit are exchanged.
- Do not disassemble the dosing unit.
- Follow the safety instructions (Safety Instructions, p. 34).

NOTICE

A functional test must always be carried out after exchanging components.

- a. Completely unscrew the retaining ring of the dispensing cartridge by hand and remove the dispensing cartridge.
- b. Insert a new dispensing cartridge into the valve block and tightly screw on the retaining ring by hand.
- c. Perform a functional test for leaks.

NOTICE**Factory adjustment**

The dispensing cartridge is pre-adjusted in the factory. Therefore, calibration is not necessary after replacement.

- Before using the instrument, rinse thoroughly and discard the first dispensed volume.
- Avoid splattering.
- Perform cleaning for trace analysis, depending on requirements (Cleaning for trace analysis, p. 50).
- Disassembly of the dispensing cartridge is not permitted!

10 Troubleshooting

Problem	Possible cause	Corrective action
Liquid is above the piston	Piston leaky	Perform standard cleaning (Standard cleaning, p. 49), replace dispensing cartridge (Replacing the dispensing cartridge, p. 52).
Piston sluggish	Crystal deposits, contamination	Stop dispensing immediately. Perform standard cleaning (Standard cleaning, p. 49).
Filling not possible	Volume setting screw at the lower stop	Set the desired volume (Dispensing, p. 42).
	Filling valve stuck	Clean the filling valve. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Dislodging a stick valve ball, p. 52). If necessary, replace the filling valve.
Dispensing not possible	Discharge valve stuck	Unscrew the discharge valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it. If necessary, replace the filling valve (Replacing the discharge valve, p. 51).

Problem	Possible cause	Corrective action
Discharge tube (with or without recirculation valve) cannot be mounted sufficiently	Discharge valve is not screwed in deeply enough	Tighten the discharge valve with the assembly tool until it meets the stop so that the threads are no longer visible.
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Screw connections loose	Tighten the valves firmly with the assembly tool
	The instrument has not been primed	Prime the instrument (Priming, p. 40).
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged	Perform cleaning (Cleaning, p. 49). Tighten the valves firmly with the assembly tool.
	Recirculation tube not connected	Install the recirculation tube (First Steps, p. 38).
Dispensed volume too low	Filling tube is loose or damaged	Perform cleaning (Cleaning, p. 49). Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Filling valve contaminated, loose or damaged	Perform cleaning (Cleaning, p. 49). Tighten the valves using the assembly tool. If necessary, replace filling valves.
Leaking liquid at the retaining ring	Dispensing cartridge loose or piston seal damaged	Tighten the retaining ring. If necessary, replace the dispensing cartridge
Leaking liquid between instrument and bottle	Recirculation tube not connected	Install the recirculation tube (First Steps, p. 38).
	Volatile reagent dispensed without sealing ring	Install the seal (Accessories, p. 43)

11 Product markings

Symbol or number	Meaning
	General warning symbol
	Refer to the operating manual
	Use eye protection
	Use hand protection

Symbol or number	Meaning
	Wear protective clothing
XXZXXXXX	Serial number
DE-M 18	The device is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the last two digits of the year in which the marking was affixed (here: 2018).
www.brand.de/ip	Patent information

12 Ordering Information

Dispensette[®] S Trace Analysis, Analog-adjustable

Volume ml	Valve spring	Without recirculation valve Order No.	With recirculation valve Order No.
10	Platinum-iridium	4640 040	4640 041
10	Tantalum	4640 240	4640 241

See scope of delivery, Scope of delivery, p. 33

13 Accessories/spare parts

Bottle adapter

External thread	for bottle thread/ground joint	Material	Order No.
GL 32	GL 25	ETFE	7043 75
GL 32	GL 28/ S 28	ETFE	7043 78
GL 32	GL 45	ETFE	7043 95
GL 45	GL 32	ETFE	7043 98
GL 45	GL 38	ETFE	7043 99
GL 45	S* 40	PTFE	7043 91

* buttress thread

Discharge tubes

With and without recirculation valve.

Nominal volume 10 ml.

Screw cap ETFE.

Discharge tube marked with 'Pt-Ir' or 'Ta'.

Packaging unit 1 piece.

Valve spring	Length mm	Without recirculation valve Order No.	With recirculation valve Order No.
Platinum-iridium	105	7080 22	7081 22
Tantalum	105	7080 24	7081 24

Flexible discharge tube with recirculation valve

PTFE, spiraled, approx. 800 mm long, with safety handle.

Packaged unit 1 pc

Not suitable for hydrofluoric acid (HF).

Nominal volume ml	Discharge tube Outer diameter	Discharge tube Inner diameter	Order No.
10	3	2	7081 3 2

Dispensette® S Trace Analysis discharge valve

PFA/sapphire

Valve marked with 'Pt-Ir' or 'Ta'

Packaged unit 1 pc

for nominal volume ml	Valve spring	Order No.
10	Platinum-iridium	6732
10	Tantalum	6733

Dispensette® S Trace Analysis filling valve

PFA/sapphire


Valve marked with 'S'

Packaged unit 1 pc

for nominal volume ml	Order No.
10	6739

Telescopic filling tube

FEP. Individually adjustable lengths. Packaged unit 1 pc

For nominal volume ml	Outer diameter mm	Length mm	Order no.
10	6	70-140	7082 10
"	"	125-240	7082 12
"	"	195-350	7082 14
"	"	250-480	7082 16

Dispensing cartridge


Dispensing cartridge with retaining ring. Nominal volume 10 ml, adjusted, includes performance certificate.

Packaged unit 1 pc

Description	Order no.
Dispensing cartridge	7080 35

Air vent cap for micro filter with Luer-cone	Recirculation tube	Adjustment/assembly tool
PP. Air vent cap and PTFE-sealing ring. Packaged unit 1 pc Order No. 7044 95	FEP Packaged unit 1 pc Order No. 6747	Packaged unit 1 pc Order No. 6748

Bottle stand	Sealing ring for valve block	Drying tube
 <p data-bbox="92 422 389 520">PP. Support rod, 325 mm, base plate 220 x 160 mm. Packaged unit 1 pc Order No. 7042 75</p>	 <p data-bbox="395 287 692 359">PTFE, for highly volatile media. Packaged unit 1 pc Order No. 7044 86</p>	 <p data-bbox="697 287 991 391">Incl. PTFE sealing ring (without drying agent). Packaged unit 1 pc Order No. 7079 30</p>

Screw cap with fastener

Packaged unit 1 pc

Description	Nominal volume ml	Order No.
ETFE	10	7060 29

14 Repairs

14.1 Sending for repair

NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech[®] Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1-860-767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, ‘C’ Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Guangqi Culture Plaza
Room 506, Building B
No. 2899, Xietu Road
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

15 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

16 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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BRAND Scientific Equipment Pvt.Ltd.