Gebrauchsanleitung | Operating manual |  
Mode d'emploi | Instrucciones de manejo |  
Istruzione | 操作手册

## HandyStep<sup>®</sup> S

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10.1 Accuracy table .....	46
10.2 Dispenser tips from other manufacturers.....	46
<b>11 Ordering Information.....</b>	<b>47</b>
<b>12 Repairs .....</b>	<b>48</b>
12.1 Sending for repair .....	48
<b>13 Calibration service.....</b>	<b>51</b>
<b>14 Warranty .....</b>	<b>52</b>
<b>15 Disposal.....</b>	<b>52</b>

# 1 Scope of delivery

HandyStep® S, DE-M marking, with serial number, quality certificate, shelf mount, 3 PD-Tips: 0.1 ml, 1 ml and 10 ml

## 2 Terms of use

- Please carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: [www.brand.de](http://www.brand.de).

### 2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

### 2.2 Symbols

Symbol	Meaning
	Hazardous area

## 2.3 Format

Layout	Meaning	Layout	Meaning
<b>1. Task</b>	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

# 3 Safety Instructions

## 3.1 General safety instructions

### Please read carefully!

The HandyStep® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual before operation.
2. Follow the general hazard instructions and safety regulations (e.g., wear protective clothing, eye protection, and protective gloves). When working with infectious or hazardous samples, the standard laboratory rules and precautions must be adhered to.
3. Observe all specifications provided by reagent manufacturers.
4. Use the instrument only for dispensing liquids within the defined limitations and restrictions of use. Comply with the operating exclusions; see Exclusions of use, p. 33. If in doubt, contact the manufacturer or supplier.
5. Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels.
6. Avoid touching the tip opening when working with aggressive media.

7. Never use force on the instrument.
8. Do not attempt to make any technical alterations. Disassembly of the instrument is not permitted.
9. Always check that the instrument is in proper working condition before use. If an instrument fault is signaled, stop dosing immediately and follow the instructions in the chapter Troubleshooting, p. 44 . Contact the manufacturer if needed.

## 3.2 Function

HandyStep® S is a repetitive pipette for fast and simple repetitive liquid dispensing. The combination of PD-Tip size and the stroke setting made on the stroke setting dial enables volumes between 2 µl and 5 ml to be dispensed with maximum precision and accuracy.

## 3.3 Usage limits

See Limitations of use, p. 34 .

## 3.4 Exclusions of use

- When the instrument is operated properly, the dosed sample comes in contact only with the tip and not with the HandyStep® S itself.
- The user is responsible for checking the suitability of the instrument for the intended use. Before using third-party dispenser tips, the user must verify their suitability and carry out a function/volume test!
- Do not use the instrument for dispensing liquids that corrode polypropylene (including cylinder PD-Tips), polyethylene (piston PD-Tips), LCP (piston 0.1 ml PD-Tip) or PC/PBT and polycarbonate.

### NOTICE


Instrument and tips are not autoclavable (sterile PD-Tips see Ordering Information, p. 47 ). PD-Tips are disposable products. Standard disinfectants can be used to disinfect the instrument, as long as they do not corrode PC/PBT and PC.

## 3.5 Limitations of use

<b>Operating range *)</b>	15 °C ... 40 °C (59 °F ... 104 °F)
<b>Vapor pressure</b>	Up to 500 mbar
<b>Viscosity</b>	20 mPa s at 50 ml PD-Tip 260 mPa s at 5 ml PD-Tip 977 mPa s at 1.25 ml PD-Tip

\*) Additional temperatures upon request

## 4 Functional and operating elements

1	Volume table (on back side of instrument under transparent clip)	5	Ejector button (release cylinder, eject PD-Tip)
2	Stroke setting dial	6	Locking/filling lever - Lock/decouple piston - Aspirate - Empty completely
3	Finger rest	7	Dispensing lever in start position
4	PD-Tip II	8	Stroke indicator (number = mm per stroke)

## Volume table

		Tip size (ml)					
Scaling		0.1	0.5	1	2.5	Steps	
1	2	10	20	50	49		
1	3	15	30	75	32		
2	4	20	40	100	24		
3	5	25	50	125	18		
3	6	30	60	150	15		
4	7	35	70	175	13		
4	8	40	80	200	11		
5	9	45	90	225	10		
5	10	50	100	250	9		
		Volume (µl)					

The volume table on back of the instrument is two-sided. To change the table, push the clip up, turn the table around and push the clip back down again.





## Shelf/rack mount

The shelf mount is comes with adhesive strips for installation. It can be used without the strip in the Transferpette® S 6-position table stand.

## 5 Operation

### 5.1 Inserting a PD-Tip

#### NOTICE

PD-Tips are disposable products!

- a. Push the locking/filling lever to the lower stop.
- b. Swing out the locking/filling lever.
- c. Insert PD-Tip straight and upright from the bottom until it clicks into place.
- d. Swing back the locking/filling lever.

## 5.2 Setting the volume

The possible dispensing volumes are listed in tables. These are found under the transparent clip on the back of the instrument.

9 half-stroke positions can be set, from 1 to 5.

- Search for the desired dispensing volume (volume) in the table.
- Use the stroke setting dial to set the stroke number for the volume that corresponds to the PD-Tip (tip size).
- Dispensing steps: depending on the PD-Tip used, the same volume can be dispensed in a varying number of steps. The smaller the number of steps, the higher the accuracy of the dispensed volume.

Example:

Desired dispensing volume per step: **200 µl**

PD-Tip size used:

**2.5 ml** = stroke setting **4** = max. **11** dispensing steps

**5 ml** = stroke setting **2** = max. **24** dispensing steps


**10 ml** = stroke setting **1** = max. **49** dispensing steps

Setting	Tip size (ml)										
	0.1	0.5	1	1.25	2.5	5	10	12.5	25	50	
1	2	10	20	25	50	100	200	250	500	1000	49
1.5	3	15	30	37.5	75	150	300	375	750	1500	32
2	4	20	40	50	100	200	400	500	1000	2000	24
2.5	5	25	50	62.5	125	250	500	625	1250	2500	19
3	6	30	60	75	150	300	600	750	1500	3000	15
3.5	7	35	70	87.5	175	350	700	875	1750	3500	13
4	8	40	80	100	200	400	800	1000	2000	4000	11
4.5	9	45	90	112.5	225	450	900	1125	2250	4500	10
5	10	50	100	125	250	500	1000	1250	2500	5000	9

Volume (µl)

**NOTICE**

The volumes listed in the table are the only options possible. Since the first dispensing step must be discarded, the number listed in the table is always 1 step less.

**5.3 Filling a PD-Tip****a.****b.****c.****NOTICE**


We recommend priming a new PD-Tip before use. To do this, fill the tip with a minimal amount of liquid and empty it again. Small air bubbles in the area of the piston after priming do not affect the results.

- a.** Push the locking/filling lever to the lower stop.
- b.** Immerse the head of the PD-Tip vertically into the liquid 3 - 10 mm.
- c.** Slowly pull the locking/filling lever up to the upper stop to prevent the formation of air bubbles.

**NOTICE**

Small air bubble near the piston can be ignored since a locking mechanism prevents the residual volume from being dispensed inadvertently after the final dispensing step.

## 5.4 Repetitive dispensing

### NOTICE

**The first dispensing step must be discarded!**

- Check the volume setting once again. (Compare with the settings in the volume table!)
- Wipe off any liquid sticking to the head of the PD-Tip with a lint-free pulp cloth.
- Place the head of the PD-Tip against the vessel wall.
- Dispense the liquid by pressing down the dispensing lever completely, then let it slide back up again.
- Make sure to dispense smoothly and evenly.

## 5.5 Ejecting a PD-Tip

### ⚠ WARNING


Residual media may be present at the head of the PD-Tip. Only eject a PD-Tip in a manner that does not endanger the user or other people.

### NOTICE

If viscous media is repeatedly dispensed, the leak tightness of the PD-Tip can no longer be guaranteed.

- a. Hold the HandyStep® S over a vessel.
  - b. Empty the PD-Tip by pushing the locking/filling lever to its lower stop.
  - c. Swing out the locking/filling lever (this will cause the piston to decouple).
  - d. Hold the HandyStep® S over a waste container. Press the ejector button down (releases the cylinder).
- ✓ PD-Tip is ejected.

## 5.6 Adapter for 25 ml and 50 ml

Adapter

To use the 25 ml and 50 ml PD-Tips, the reusable adapter supplied with the tips is required. This is coupled to the PD-Tips with a bayonet lock and can be detached after use. Adapters are autoclavable at 121 °C (2 bar), in accordance with DIN EN 285.

## 6 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to correspond with individual requirements. Gravimetric volume testing according to DIN EN ISO 8655-5 is performed as follows:

### 1. Preparing the instrument

Insert a tip. This test can be carried out with PD-Tips of any size. However, a 5 ml PD-Tip is normally used for this purpose.

### 2. Performing the test

- a.** Adjust the stroke setting of the HandyStep®S to Step 5 (10% of nominal volume).
- b.** Fill the PD-Tip. Immerse the PD-Tip vertically into the test liquid.
- c.** The first dispensing step must be discarded. It serves to compensate for play.
- d.** Place a weighing vessel filled with deionized water on the scale and tare the scale.
- e.** Dispense the second step into the weighing vessel. In doing so, push the dispensing lever down at a uniform speed until the first stop and hold it. Then wipe the tip off over a distance of approx. 10 mm.
- f.** Enter the weight into the test record.
- g.** Repeat points D to F a total of 10 times. (For the 10th volume test, the 5 ml PD-Tip with a stroke setting of 5 must be refilled – points A to C).

**h.** Repeat the same testing procedure at stroke setting 3 (6% of nominal volume) and 1 (2% of nominal volume).

**i.** This results in a total of the 30 required weighed values.

$x_i$  = Weighing results                       $n$  = Number of weighings

$Z$  = Correction factor  
(e.g. 1.0029  $\mu\text{l}/\text{mg}$  at 20 °C, 1013 hPa)

$$\text{Mean} \quad \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume} \quad \bar{V} = \bar{x} \cdot z$$

$$\text{Accuracy*} \quad A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

$V_0$  = Nominal volume

$$\text{Coefficient of variation* CV\%} = \frac{100 s}{\bar{V}}$$

$$\text{Standard deviation*} \quad s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

\*) = Calculation for accuracy (A%) coefficient of variation (CV%):

A% and CV% are calculated using the formulas of statistical quality control.

### NOTICE

Test instructions (SOPs) and a version of the calibration software EASYCAL™ 4.0 are available for download at [www.brand.de](http://www.brand.de).

## 7 Cleaning and disinfection

### 7.1 Cleaning

#### NOTICE

#### The instrument is not autoclavable

The instrument is factory-calibrated and maintenance-free.

When the instrument is dirty, clean the outer surface with a moist cloth. Use water or a mild detergent solution. Avoid acidic or aggressive cleaners.

Disassembly of the instrument is not permitted.

### 7.2 UV disinfection

The device is resistant to normal exposure to a UV disinfection lamp. The device may change color due to the UV exposure.

UV disinfection recommendations:

Light spectrum	UV-C
Wave length	220 nm ... 270 nm
Exposure time per cm <sup>2</sup> and intended inactivation rate	2 s ... 300 s

## 8 Troubleshooting

Problem	Possible cause	Corrective action
PD-Tip piston cannot be locked after the PD-Tip has been inserted.	Filling/locking lever is not pushed down to bottom position and is not swung out.	First press the ejector button and remove the PD-Tip, then push the filling/locking lever all the way down and swing it out completely.
Undefined dispensing volume	The stroke setting dial is not properly engaged.	Ensure that the stroke setting dial is securely locked into place.

Problem	Possible cause	Corrective action
Air bubble under the PD-Tip piston increasing in size.	PD-Tip is leaking.	Replace PD-Tip.
A and/or CV exceed tolerances.	PD-Tip is leaking.	Replace PD-Tip.
	Dispenser tips are difficult to move.	Replace PD-Tip or use other more pliant dispenser tips.
	Dispensing lever has not been pressed completely and with even pressure.	Press the dispensing lever with even pressure to the stop position.
	The first dispensing step was not discarded.	Discard the first dispensing step.
PD-Tip cannot be decoupled after emptying completely.	The PD-Tip has been inserted incorrectly.	Press the ejector button and then push the filling/locking lever all the way down and swing it out.

## 9 Product markings

Symbol or number	Meaning
	General warning symbol
	Follow the operating instructions
XXZXXXXX	Serial number
	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed (here: 2018).

## 10 Technical data

### 10.1 Accuracy table

PD-Tips // from BRAND, 20 °C 'Ex'

PD-Tips sizes	Volume range (µl)	A* ± % stroke setting % of nominal volume			CV* ≤ % stroke setting % of nominal volume		
		1 2%	3 6%	5 10%	1 2%	3 6%	5 10%
0.1 ml	2 - 10	4.0	2.4	1.6	6.0	3.0	2.0
0.5 ml	10 - 50	2.5	1.5	1.0	2.5	1.5	1.0
1.0 ml	20 - 100	2.5	1.5	1.0	2.0	1.2	0.8
1.25 ml	25 - 125	2.5	1.4	0.9	2.0	1.1	0.7
2.5 ml	50 - 250	1.8	1.1	0.7	1.5	0.9	0.6
5.0 ml	100 - 500	1.8	1.1	0.7	1.5	0.9	0.7
10.0 ml	200 - 1000	1.8	1.1	0.7	2.0	1.2	0.8
12.5 ml	250 - 1250	1.8	1.1	0.8	3.2	2.0	1.4
25.0 ml	500 - 2500	1.5	0.9	0.6	3.0	1.5	1.0
50.0 ml	1000 - 5000	1.5	0.8	0.5	5.0	1.8	1.2

A\* = Accuracy, CV = Coefficient of Variation

Error limits in relation to the set dispensing volume depending on the PD-Tip size, at equal temperature (20 °C) of instrument, tip, ambient environment and distilled water, and with uniform, jolt-free handling. The test is carried out in accordance with DIN EN ISO 8655-5. The nominal volume is the maximum volume printed on the PD-Tip.

### 10.2 Dispenser tips from other manufacturers

We recommend using BRAND PD-Tips //; however, the HandyStep® S can also be used with compatible dispenser tips from other manufacturers. When using compatible dispenser tips, the user must verify their suitability and carry out a function/volume test!

# 11 Ordering Information

## HandyStep S

Description	Packaging unit	Order no.
HandyStep® S, 1 shelf mount, 3 PD-Tips II (0.1 ml, 1 ml, 10 ml),	1	7051 10

## Shelf/rack mount

Description	Packaging unit	Order no.
Shelf/rack mount	1	7051 30

## PD-Tip II

Volume [ml]	Order No.	Packing unit [pcs.]	Order No. BIO-CERT® LIQUID HANDLING STERILE	Packing unit [pcs.]
0.1	7057 00	100	7057 30	100
0.5	7057 02	100	7057 32	100
1	7057 04	100	7057 34	100
1.25	7057 06	100	7057 36	100
2.5	7057 08	100	7057 38	100
5	7057 10	100	7057 40	100
10	7057 12	100	7057 42	100
12.5	7057 14	100	7057 44	100
25	7057 16	25+1 Adapter	7057 46	25+1 Adapter
50	7057 18	50+1 Adapter	7057 48	25+1 Adapter
Set PD-Tips II 0,5 ml ... 12.5 ml	7057 20	per 20	—	—

## PD-Tip II Set

Description	Order no.
PD-Tip II Set, 20 PD-Tips II each in sizes 0.5, 1.0, 1.25, 2.5, 5.0, 10 and 12.5 ml	7057 20

## Adapter for 25 ml and 50 ml PD-Tip II

Volume [ml]	Order No.	Packaging unit	Property
25 ml and 50 ml	7023 98	10	
25 ml and 50 ml	7023 99	5	BIO-CERT® LIQUID HANDLING STERILE

# 12 Repairs

## 12.1 Sending for repair

### NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

### Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

## Within Germany

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

## Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

## Outside Germany, the U.S. and Canada

Please contact a service partner, dealer or branch office in your region. For a list of service partners and dealers, please visit [www.brand.de](http://www.brand.de). If you live in a region without an official BRAND partner, please contact our main factory in Germany **before** you sent the instrument.

## Contact addresses

### Germany:

BRAND GMBH + CO KG  
Otto-Schott-Str. 25  
97877 Wertheim (Germany)  
T +49 9342 808 0  
F +49 9342 808 98000  
info@brand.de  
www.brand.de

### India:

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park,  
Powai  
Mumbai-400 076 (India)  
T +91 22 42957790  
F +91 22 42957791  
info@brand.co.in  
www.brand.co.in

### USA and Canada:

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)  
T +1-860-767 2562  
F +1-860-767 2563  
info@brandtech.com  
www.brandtech.com

### China:

BRAND (Shanghai) Trading Co., Ltd.  
Guangqi Culture Plaza  
Room 506, Building B  
No. 2899, Xietu Road  
Shanghai 200030 (P.R. China)  
T +86 21 6422 2318  
F +86 21 6422 2268  
info@brand.com.cn  
www.brand.cn.com

## 13 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The detailed testing instructions are available for download on [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or by the BRAND DAkkS laboratory. Simply send us the instruments to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days together with a test report (factory calibration) or a DAkkS calibration certificate. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on [www.brand.de](http://www.brand.de) (see Technical Documents).

### **For customers outside Germany**

If you would like to use our calibration service, please contact one of our service partners in your region. If a factory calibration is requested, the service partner can forward the instrument to BRAND.

## 14 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

### **USA and Canada:**

Find more warranty information on [www.brandtech.com](http://www.brandtech.com).

## 15 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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