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# Table of contents

|          |  |           |
|----------|--|-----------|
| <b>1</b> | <b>Introduction .....</b>                            | <b>47</b> |
| 1.1      | Scope of delivery.....                               | 47        |
| 1.2      | Use the instructions for use.....                    | 47        |
| <b>2</b> | <b>Safety regulations.....</b>                       | <b>48</b> |
| 2.1      | Safety regulations .....                             | 48        |
| 2.2      | Target group.....                                    | 49        |
| 2.3      | Purpose .....  | 49        |
| 2.4      | Use .....  | 50        |
| 2.5      | Improper use .....                                   | 50        |
| 2.6      | Foreseeable misuse .....                             | 50        |
| 2.7      | Operating Limitations.....                           | 50        |
| 2.8      | Limitations of use.....                              | 50        |
| 2.9      | Operating exclusions .....                           | 51        |
| <b>3</b> | <b>Functional and operational components .....</b>   | <b>52</b> |
| <b>4</b> | <b>Pipetting .....</b>                               | <b>56</b> |
| <b>5</b> | <b>Checking the volume .....</b>                     | <b>60</b> |
| <b>6</b> | <b>Accuracy table .....</b>                          | <b>62</b> |
| <b>7</b> | <b>Adjusting – User and factory adjustment .....</b> | <b>64</b> |
| 7.1      | Permanent factory adjustment: Easy Calibration ..... | 64        |
| 7.2      | Temporary adjustment: User Adjustment .....          | 66        |
| <b>8</b> | <b>Disinfection/autoclaving.....</b>                 | <b>70</b> |
| 8.1      | Autoclaving.....                                     | 70        |
| 8.2      | UV sterilization .....                               | 70        |
| 8.3      | PE filter .....                                      | 70        |
| <b>9</b> | <b>Maintenance.....</b>                              | <b>72</b> |
| 9.1      | Cleaning – volumes up to 1,000 µL.....               | 73        |
| 9.2      | Cleaning – volumes of 2.5, 5, and 10 ml .....        | 74        |

|   |           |
|---|-----------|
| <b>10 Troubleshooting .....</b>                             | <b>75</b> |
| <b>11 Product markings .....</b>                            | <b>76</b> |
| <b>12 Order Information.....</b>                            | <b>77</b> |
| 12.1    Laboratory instrument and accessories .....         | 77        |
| 12.2    Spares .....  | 78        |
| 12.3    Additional accessories .....                        | 79        |
| <b>13 Repairs.....</b>                                      | <b>81</b> |
| 13.1    Sending for repair .....                            | 81        |
| <b>14 Calibration service .....</b>                         | <b>83</b> |
| <b>15 Information about your laboratory instrument.....</b> | <b>84</b> |
| <b>16 Warranty .....</b>                                    | <b>85</b> |
| <b>17 Disposal .....</b>                                    | <b>86</b> |

# 1 Introduction

## 1.1 Scope of delivery

Transferpette® pro, variable type, marked DE-M, with quality certificate, shelf mount, and silicone grease.

## 1.2 Use the instructions for use

- Carefully read the operating manual before using the device for the first time.
- Keep the user manual easily accessible. It is part of the device.
- Be sure to include the operating manual if you transfer possession of this device to a third party.

### 1.2.1 Signal words and their meaning

| Signal word                                 | Meaning   |
|---|---|
| <b>⚠ WARNING</b><br>or <b>⚠WARNING!</b> ... | WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.      |
| <b>⚠ CAUTION</b><br>or <b>⚠CAUTION!</b> ... | CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.     |
| <b>NOTICE</b><br>or <b>NOTICE</b> ...       | NOTE is used to address actions that are not related to physical injury. Example: Possible property damage. |

### 1.2.2 Presentation of descriptions of actions

| Format         | Meaning                                   |
|----------------|---|
| <b>1. Task</b> | Indicates a task.                         |
| a., b., c.     | Indicates the individual steps of a task. |
| >              | Indicates a prerequisite for a task.      |
| ⇒              | Indicates a result of a completed task.   |

## 2 Safety regulations

### 2.1 Safety regulations

#### Please read carefully!

The laboratory device Transferpette® pro can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the user manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/media must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.
- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the "Troubleshooting" section in the user manual. Contact the manufacturer, if necessary.

- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the user manual!
- Only autoclave the laboratory device when this is permissible according to the user manual.

## 2.2 Target group

The operating manual is intended for users who use the laboratory instrument in the course of their professional activities. Users are familiar with the typical safety regulations and working methods in laboratories and have been trained accordingly. They can recognize potential hazards and protect themselves from them. The operating manual assumes this expertise and does not replace basic laboratory training or specific safety training.

## 2.3 Purpose

This is an air displacement pipette for pipetting liquids of medium density and low to medium viscosity.

## 2.4 Use

Use the laboratory device Transferpette® pro only for pipetting or dispensing liquids within the defined limits of use.

## 2.5 Improper use

Various risks may arise if the laboratory instrument is used improperly. These risks include: inaccurate liquid delivery, damage to the laboratory instrument, and the risk of contamination, infection, and injury from contact with the pipetted media.

Any use other than for pipetting or dispensing liquids within the defined operating limits is considered improper use.

## 2.6 Foreseeable misuse

A typical misuse is pipetting or dispensing liquids with too high viscosity or using unsuitable tips.

## 2.7 Operating Limitations

Volume accuracy may be affected by liquid density, special tip shapes, or deviations from ambient temperature.

However, the temporary user adjustment enables you to correct any resulting deviations and improves accuracy under conditions that differ from the factory setting (aqueous medium, ISO 8655). See Temporary adjustment: User Adjustment, p. 66 .

## 2.8 Limitations of use

The pipette is used for dispensing liquids within the following limits:

- Operating temperature from +15 °C to +40 °C (59 °F to 104 °F). Additional temperatures upon request
- vapor pressure up to 500 mbar
- Viscosity: 260 mPa s

With viscous liquids, the dispensing speed may need to be adjusted.

## 2.9 Operating exclusions


The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:

- for liquids with very high steam pressure
- Liquids that corrode the following materials:
  - Fluoroelastomer rubber (FKM)
  - Polyamide (PA)
  - Polycarbonate (inspection window)
  - Polyetheretherketone (PEEK)
  - Polyphenylene sulfide (PPS)
  - Polypropylene (PP)
  - Polyvinylidene fluoride (PVDF)

Additional information on the chemical resistance of plastics can be found at [www.brand.de](http://www.brand.de).

## 3 Functional and operational components

### Front

|                              |  |
|------------------------------|--|
| Adjustable finger rest       | The laboratory instrument Transferpette® pro has an adjustable finger rest. This allows you to adapt the pipette to your preferred grip, see " Pipetting".<br><br>The instrument can be labeled at the finger rest: to do so, remove the label window from the finger rest and take out the label strip.       |
| Volume-change protection     | The volume adjustment lock prevents the volume from being changed during pipetting work. To unlock it, slide the volume adjustment lock toward the pipetting button.   |
| Volume display               | The numbers in the display are read from top to bottom; the white dash corresponds to the decimal point.   |
| Counter mechanism protection | Once the volume adjustment lock is released, set the volume using the volume-setting wheel. The integrated counter mechanism protection overrides the volume-setting wheel once the maximum or minimum volume setting is reached: the volume-setting wheel remains rotatable but no longer adjusts the volume. |
| Hand grip                    | Screw the pipetting unit into the hand grip. Attach the tip onto the tip cone.   |

## Rear

Permanent adjustment to factory settings:  
Easy Calibration

Temporary adjustment to  
changing conditions:  
User Adjustment

User Adjustment scale

Cover

Seal

Label window

Serial number and  
product markings

**QR code:**  
Use your smartphone, tablet,  
or webcam to scan and visit  
the following website:  
[www.brand.de/myproduct](http://www.brand.de/myproduct)

The website contains  
serialized information about  
your laboratory instrument.

## Easy Calibration Technology

The Easy Calibration function is located in the finger rest and is used to reset the pipette to its factory settings (see Adjusting – User and factory adjustment, p. 64 ).

## User Adjustment Technology

The hand grip also includes the User Adjustment function. This allows the pipette to be adjusted for specific liquids and dispensing conditions. See Temporary adjustment: User Adjustment, p. 66

The switch for user adjustment is located behind the cover. A seal is applied at delivery. Remove it upon first use and dispose of it.

## QR Code and serialized information

The QR code provides access to [www.brand.de/myproduct](http://www.brand.de/myproduct) and serialized information about your pipette.

If you wish to access the MyProduct information without the QR code, you will also need the order number and serial number of your pipette.

## 4 Pipetting

a.

- a. Turn the finger rest to a comfortable working position.

b.

- b. Slide the volume adjustment lock in the direction shown, against slight resistance.

c.

- c. Set the volume using the volume-setting wheel.

**NOTICE** If the volume-setting wheel is turned beyond the maximum or minimum volume, it glides over the volume adjustment and thus protects the counter mechanism from damage.

d.

- d. Close the volume adjustment lock.

⇒ The volume-setting wheel can still be turned but will no longer change the volume.

e.

- e. Attach the pipette tip. Ensure a secure fit.

The 2 ml to 10 ml pipettes should only be used with a built-in PE filter (see UV sterilization, p. 70 ).

**NOTICE** Pipette tips are disposable products!


f.

- f. Press the pipetting button down to the first stop.

**NOTICE** We recommend rinsing the pipette tip five times with the liquid (aspirate and dispense) before pipetting in order to ensure maximum precision and accuracy.

g.

g. Immerse the tip in the liquid.

h.

h. Slowly release the pipetting button.

⇒ Liquid is aspirated.

Keep the tip submerged until the volume is fully aspirated. Increase the wait time when pipetting more viscous liquids or larger volumes.

**⚠ CAUTION!** Do not lay the pipette down with filled tips. Contamination may occur!

| Volume range | Immersion depth [mm] | Wait time [s] |
|--------------|----------------------|---------------|
| 0.1–1 µl     | 1–2                  | 1             |
| 1–100 µl     | 2–3                  | 1             |
| 100–1000 µl  | 2–4                  | 1             |
| > 1,000 µl   | 3–6                  | 3             |

i.

- i. To dispense the liquid, hold the tip at an angle against the vessel wall, slowly press the pipetting button, and wipe the tip.

To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density fluids.

To fully empty the tip, press the pipetting button down to the second stop (Fig. f).

j.

- j. To remove the tip, hold the pipette over a container and press the tip ejection button.

## Storage

You can also hook the Transferpette® pro in the holder or rack with an adjustable finger rest.

**⚠ CAUTION!** Do not hang the pipette with a filled tip in the holder.  
Contamination may occur!

## 5 Checking the volume

We recommend testing the instrument every three to 12 months depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. The complete testing procedure (SOP) can be downloaded at [www.brand.de](http://www.brand.de).

You can download the detailed test instructions (SOP) at [www.brand.de](http://www.brand.de). For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from <https://shop.brand.de/>.

Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655:2022.

- a. Set User Adjustment to 0 (for procedure, see Temporary adjustment: User Adjustment, p. 66 )
- b. Set the maximum specified instrument volume (for procedure, see Pipetting, p. 56 ).
- c. Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.
- d. Aspirate the test liquid and pipette into the weighing vessel.
- e. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- f. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- g. At least 10 pipetting series and weighings in three volume ranges (100%, 50%, 10%) are recommended. Two tips must be used for each volume range to be tested.

### Calculation (for nominal volume)

$x_i$  = weighing results

$n$  = number of weighings

$V_0$  = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

**Mean:**

$$\bar{x} = \frac{\sum x_i}{n}$$

**Mean volume:**

$$\bar{V} = \bar{x} * Z$$

**Accuracy\*:**

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

**Standard deviation\*:**

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

**Coefficient of variation\*:**

$$CV\% = \frac{100 s}{\bar{V}}$$

\*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

## 6 Accuracy table

| Volume range [µl] | Partial volume<br>[µl] | A ≤ ± % | CV ≤ ± % | Tip type [µl] |
|-------------------|------------------------|---------|----------|---------------|
| 0.1-1             | 1                      | 2       | 1.2      | 0.1-20        |
|                   | 0.5                    | 4       | 2.4      |               |
|                   | 0.1                    | 20      | 12       |               |
| 0.1-2.5           | 2.5                    | 1.4     | 0.7      | 0.5-20        |
|                   | 1.25                   | 2.5     | 1.5      |               |
|                   | 0.25                   | 12      | 6        |               |
| 0.5-10            | 10                     | 1       | 0.5      | 0.5-20        |
|                   | 5                      | 1.6     | 1        |               |
|                   | 1                      | 7       | 4        |               |
| 2-20 gray         | 20                     | 0.8     | 0.4      | 0.5-20        |
|                   | 10                     | 1.2     | 0.7      |               |
|                   | 2                      | 5       | 2        |               |
| 2-20 yellow       | 20                     | 0.8     | 0.4      | 2-200         |
|                   | 10                     | 1.2     | 0.7      |               |
|                   | 2                      | 5       | 2        |               |
| 5-50              | 50                     | 0.8     | 0.3      | 2-200         |
|                   | 25                     | 1.2     | 0.5      |               |
|                   | 5                      | 4       | 2        |               |
| 10-100            | 100                    | 0.6     | 0.2      | 2-200         |
|                   | 50                     | 0.8     | 0.4      |               |
|                   | 10                     | 3       | 1        |               |
| 20-200            | 200                    | 0.6     | 0.2      | 2-200         |
|                   | 100                    | 0.8     | 0.3      |               |
|                   | 20                     | 3       | 0.6      |               |
| 30-300            | 300                    | 0.6     | 0.2      | 5-300         |
|                   | 150                    | 0.8     | 0.3      |               |
|                   | 30                     | 3       | 0.6      |               |

| Volume range [µl] | Partial volume<br>[µl] | A ≤ ± % | CV ≤ ± % | Tip type [µl] |
|-------------------|------------------------|---------|----------|---------------|
| 100-1000          | 1000                   | 0.6     | 0.2      | 50-1000       |
|                   | 500                    | 0.8     | 0.3      |               |
|                   | 100                    | 3       | 0.6      |               |
| 250-2500          | 2500                   | 0.6     | 0.2      | 500-5000      |
|                   | 1250                   | 0.8     | 0.3      |               |
|                   | 250                    | 3       | 0.6      |               |
| 500-5000          | 5000                   | 0.6     | 0.2      | 500-5000      |
|                   | 2500                   | 0.8     | 0.3      |               |
|                   | 500                    | 3       | 0.6      |               |
| 1000-10000        | 10000                  | 0.6     | 0.2      | 1000-10000    |
|                   | 5000                   | 0.8     | 0.3      |               |
|                   | 1000                   | 3       | 0.6      |               |

\*A = Accuracy, CV = Coefficient of Variation

Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

## 7 Adjusting – User and factory adjustment

You have the following options to adjust the instrument:

- Factory adjustment:

The factory adjustment is used for permanent calibration of the instruments to aqueous media according to ISO 8566 in cases of volume deviations.

- Temporary User Adjustment:

The User Adjustment is used for temporary volume adjustment under changing conditions. It can be reset to the original state (factory adjustment).

### 7.1 Permanent factory adjustment: Easy Calibration

The instrument is factory-calibrated to aqueous solutions in accordance with ISO 8655. If it is determined that the pipette is inaccurate, it can be adjusted using the Easy Calibration technique.

a.

- a. Check whether the User Adjustment is set to 0 (see Temporary adjustment: User Adjustment, p. 66 ).

**NOTICE** If the User Adjustment is not set to 0, the pipette will be misaligned when attempting to perform factory adjustment. In this case, set the User Adjustment to 0 and repeat the factory adjustment as described.

- b. Perform a volume check and determine the actual value (see Checking the volume).

**c.**

- c. Slightly lift and set aside the labeling window (1) on the finger rest with your thumb. Use a paperclip or an unused pipette tip to remove and dispose of the protective foil (2).

**d.**

- d. Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

**e.**

- e. Set the volume adjustment lock to the UNLOCK position and adjust the previously determined actual volume using the volume-setting wheel. Position UNLOCK (see Pipetting, p. 56 > “Set volume”). Set the volume adjustment lock back to the LOCK position. After each adjustment, a volume check is recommended.

f.

- f. Slide the adjustment slider completely back again, allow the volume-setting wheel to slide downward and release the adjustment slider. If the volume-setting wheel does not slide down easily, move it slightly back and forth until it clicks into place. Reinsert the label window.

**NOTICE** The change to factory settings is indicated by the red adjustment slider now visible in the label window.

## 7.2 Temporary adjustment: User Adjustment

Temporary User Adjustment improves accuracy under conditions that deviate from the factory settings (aqueous medium, ISO 8655). This enables temporary volume corrections under changing conditions because deviations from water in physical properties, temperature differences between liquid and ambient conditions, specific tip designs, and environmental factors can all affect accuracy.

**NOTICE** User Adjustment modifies the volume setting by a certain volume offset (e.g. 100 µl: +2 µl = +2%). If the volume setting changes significantly (e.g. 10 µl: +2 µl = +20%), the adjustment value must be recalculated.

## Setting the User Adjustment

- Pry off and remove the cover (1) and seal (2) (e.g., using a paperclip). Dispose of the seal.
- Slide the slider (3) down into the recess and hold it there. Use the volume-setting wheel (4) to set the desired User Adjustment value (see below) on the scale. Release the volume-setting wheel and slowly return the slider (3).

**NOTICE** If the slider is stuck, gently push it back toward the recess (3) and slowly return it again.

- ⇒ The value is set when the user adjustment value aligns with the marking (5).
- Reinsert the cover (1).
  - Verify the adjustment gravimetrically.

## Determining the User Adjustment

Example: Pipetting 180 µl with a 20–200 µl pipette

- Perform control weighings on a precision balance and calculate the actual volume:  
Actual volume: 178.4 µl
- Calculate the volume to be corrected:  
Volume offset: 1.6 µl (= 180 µl – 178.4 µl)

c. Determine and set the User Adjustment value:

For our 200 µl instrument, each dash corresponds to a step value of 0.2 µl (see assignment table). A volume offset of 1.6 µl is added by setting to +8 (= 1.6 µl / 0.2).

$$\text{Actual volume} = \frac{\text{Mean of liquid weights}}{\text{Density of liquid} - \text{Density of air (0.0012 g/ml)}}$$

$$\text{Volume offset} = \text{Target volume} - \text{Actual volume}$$

$$\text{User Adjustment value} = \frac{\text{Volume offset}}{\text{Step value}}$$

### Assignment table for User Adjustment

| Nominal volume<br>[µl] | The highlighted column [1] indicates the step value for the respective instrument. |       |        |       |        |        |   |       |       |      |       |      |       |      |       |
|------------------------|--|-------|--------|-------|--------|--------|---|-------|-------|------|-------|------|-------|------|-------|
|                        | -25  | -20   | -15    | -10   | -5     | -1     | 0 | 1     | 5     | 10   | 15    | 20   | 25    | 30   | 35    |
| 1                      | -0,025   | -0,02 | -0,015 | -0,01 | -0,005 | -0,001 | 0 | 0,001 | 0,05  | 0,01 | 0,015 | 0,02 | 0,025 | 0,03 | 0,035 |
| 2,5                    | -0,05  | -0,04 | -0,03  | -0,02 | -0,01  | -0,002 | 0 | 0,002 | 0,01  | 0,02 | 0,03  | 0,04 | 0,05  | 0,06 | 0,07  |
| 10                     | -0,25  | -0,2  | -0,15  | -0,1  | -0,05  | -0,01  | 0 | 0,01  | 0,05  | 0,1  | 0,15  | 0,2  | 0,25  | 0,3  | 0,35  |
| 20                     | -0,5   | -0,4  | 0,3    | 0,2   | -0,1   | -0,02  | 0 | 0,02  | 0,1   | 0,2  | 0,3   | 0,4  | 0,5   | 0,6  | 0,7   |
| 50                     | -1,25  | -1    | -0,75  | -0,5  | -0,25  | -0,05  | 0 | 0,05  | 0,25  | 0,5  | 0,75  | 1    | 1,25  | 1,5  | 1,75  |
| 100                    | -2,5   | -2    | -1,5   | -1    | -0,5   | -0,1   | 0 | 0,1   | 0,5   | 1    | 1,5   | 2    | 2,5   | 3    | 3,5   |
| 200                    | -5   | -4    | -3     | -2    | -1     | -0,2   | 0 | 0,2   | 1     | 2    | 3     | 4    | 5     | 6    | 7     |
| 300                    | -6,225   | -4,98 | -3,735 | -2,49 | -1,245 | -0,249 | 0 | 0,249 | 1,245 | 2,49 | 3,735 | 4,98 | 6,225 | 7,47 | 8,715 |
| 1000                   | -25  | -20   | -15    | -10   | -5     | -1     | 0 | 1     | 5     | 10   | 15    | 20   | 25    | 30   | 35    |
| 1250                   | -25  | -20   | -15    | -10   | -5     | -1     | 0 | 1     | 5     | 10   | 15    | 20   | 25    | 30   | 35    |
| 2500                   | -50  | -40   | -30    | -20   | -10    | -2     | 0 | 2     | 10    | 20   | 30    | 40   | 50    | 60   | 70    |
| 5000                   | -125   | -100  | -75    | -50   | -25    | -5     | 0 | 5     | 25    | 50   | 75    | 100  | 125   | 150  | 175   |
| 10000                  | -250   | -200  | -150   | -100  | -50    | -10    | 0 | 10    | 50    | 100  | 150   | 200  | 250   | 300  | 350   |

Volume offset for excess volume

Volume offset for missing volume

**NOTICE** The table shows the mechanical relationship between the steps of the User Adjustment. The volume changes indicated are approximate values and apply to the entire volume range of the instrument.

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### Recognizing user adjustment

If the red switch is visible on the back of the instrument, it has already been adjusted by a user with the user adjustment. Check whether this adjustment still suits your application (e.g., by performing a control weighing of the pipetted volume). Reset User Adjustment if necessary.

### Restore factory adjustment, reset User Adjustment

To reset the User Adjustment, set it to 0 on the scale. This restores the factory adjustment state. We recommend performing a volume check afterward.

## 8 Disinfection/autoclaving

### 8.1 Autoclaving

**NOTICE Carry out a self-test of the effectiveness of the autoclave.**

Maximum safety is achieved through vacuum sterilization. We recommend using sterilization bags.

- a. Eject the pipette tip.
- b. Pack the instrument in a sterilization bag, observing any applicable packaging regulations.
- c. Autoclave the complete pipette without further disassembly.  
Recommendation for autoclaving according to DIN EN 285 (see table below).
- d. Allow the pipette to completely cool and dry.

|                            |        |
|----------------------------|--------|
| Temperature                | 121 °C |
| Pressure                   | 2 bar  |
| Holding time in autoclaves | 15 min |

After autoclaving, tighten the connection between the hand grip and the pipette shaft if necessary.

### 8.2 UV sterilization

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

### 8.3 PE filter

PE filter for Transferpette® pro 2 ml, 2,5 ml, 5 ml, and 10 ml:

A hydrophobic PE filter is used to prevent liquid from entering the pipette.

Change the filter if it becomes wet or contaminated.

- a. Use a flat object, such as a screwdriver.
- b. Remove the filter carefully, without damaging the tip cone.

Remove the filter before autoclaving!


The device can also be operated without a filter.

## 9 Maintenance

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.


We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see Troubleshooting, p. 75 .

## 9.1 Cleaning – volumes up to 1,000 µL

- a. Detach the pipette shaft (S) from the hand grip by unscrewing it.
  - b. Unscrew the upper part of the ejector unit (A) from the pipette shaft.
  - c. Pull out the shaft (B, C and D) from the lower part of the ejector unit (E).
  - d. Unscrew the piston unit (B).
- NOTICE** Do not disassemble the piston unit (B) any further!
- e. Remove the seal with spring (C) (not possible on Transferpette® S 1 µL, 2.5 µL and 10 µL!).
  - f. Clean the parts shown with a soap solution or isopropanol, and then rinse with distilled water.
  - g. Dry the parts (max. 120 °C/248 °F).
  - h. Grease piston and seal with a very thin layer of supplied silicone grease.
- Reassemble the cooled parts in reverse order. Only hand-tighten the piston unit and the upper part of the ejector unit (A, B).

## 9.2 Cleaning – volumes of 2.5, 5, and 10 ml

- Remove the entire shaft (S) from the hand grip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
  - Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
  - Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
  - Remove the O-ring-seal from the piston unit and clean it.
- NOTICE** Do not disassemble the piston unit (G) any further!
- Clean the piston unit (G) and the bottom part of the shaft (H) with a soap solution or isopropanol, and then rinse with distilled water.
  - Dry the parts (max. 120 °C/248 °F) and allow them to cool.
  - Carefully lubricate the inside and outside of the O-ring (O) and mount it on the piston.

Reassemble the individual components in reverse order.

## 10 Troubleshooting

| <b>Problem</b>  | <b>Possible cause</b>   | <b>Corrective action</b>           |
|---|---|------------------------------------|
| Tip dripping (device leaking)   | Unsuitable tip  | Only use high-quality tips         |
|   | Tip not seated tightly  | Firmly press tip on                |
| The instrument does not aspirate or aspirates too little; the dispensed volume is too low | Seal contaminated   | Clean seal                         |
|   | Seal or cone is damaged   | Replace seal or shaft              |
|   | Piston is contaminated or damaged   | Clean or replace piston            |
| Aspiration is very slow   | Shaft is clogged  | Clean shaft                        |
| Pipette is misadjusted  | Calibrated with altered User Adjustment                                       | Set User Adjustment to 0. Readjust |
| Dispensed volume too large  | Pipetting button pressed too far (to the over-stroke point) before aspirating | Ensure proper handling.            |
| Piston sluggish   | Piston is contaminated or not greased   | Clean piston and apply grease      |

## 11 Product markings

| Symbol or number  | Meaning   |
|---|---|
|  | Read the user manual.   |
| XXZXXXXX  | Serial number   |
| DE-M 25   | The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation.<br>Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed. |
| (121 °C)  | Autoclavable up to the temperature shown  |
| Data matrix code or quick response code   | The codes link to the BRAND MyProduct website.  |
| <a href="http://www.brand.de/ip">www.brand.de/ip</a>                              | Hyperlink to BRAND patent site  |

## 12 Order Information

### 12.1 Laboratory instrument and accessories


| Nominal volume   | Order No.              |
|------------------|------------------------|
| 0.1–1 µl         | <a href="#">706868</a> |
| 0.1–2.5 µl       | <a href="#">706869</a> |
| 0.5–10 µl        | <a href="#">706870</a> |
| 2–20 µl (gray)   | <a href="#">706871</a> |
| 2–20 µl (yellow) | <a href="#">706872</a> |
| 5–50 µl          | <a href="#">706873</a> |
| 10–100 µl        | <a href="#">706874</a> |
| 20–200 µl        | <a href="#">706878</a> |
| 30–300 µl        | <a href="#">706879</a> |
| 100–1,000 µl     | <a href="#">706880</a> |
| 250–2500 µl      | <a href="#">706881</a> |
| 500–5,000 µl     | <a href="#">706882</a> |
| 1,000–10,000 µl  | <a href="#">706884</a> |

### Accessories

|   | Description   | Order No.              |
|---|---|------------------------|
|   | Table stand for 6 single-channel or multi-channel pipettes. | <a href="#">704807</a> |
|   | Table stand for 1 single-channel or multi-channel pipette.  | <a href="#">703440</a> |
|  | Wall mount  | <a href="#">704812</a> |
|  | Shelf/rack mount  | <a href="#">704811</a> |

## 12.2 Spares

### 12.2.1 Spare parts – volumes up to 1,000 µl


| Vol.             | A      | B      | C      | D       | E      |
|------------------|--------|--------|--------|---------|--------|
| 0.1–1 µl         | 705513 | 704600 | —      | 704718* | 704705 |
| 0.1–2.5 µl       | 705513 | 704667 | —      | 704717  | 704706 |
| 0.5–10 µl        | 705513 | 704601 | —      | 704721* | 704707 |
| 2–20 µl (gray)   | 705513 | 704602 | 704610 | 704727  | 704710 |
| 2–20 µl (yellow) | 705513 | 704602 | 704610 | 704723  | 704710 |
| 5–50 µl          | 705513 | 704615 | 704617 | 704722  | 704711 |
| 10–100 µl        | 705513 | 704654 | 704661 | 704724  | 704712 |
| 20–200 µl        | 705513 | 704655 | 704662 | 704725  | 704713 |
| 30–300 µl        | 705513 | 704668 | 704664 | —       | 704714 |
| 100–1000 µl      | 705513 | 704656 | 704663 | 704726  | 704715 |

\* Seal permanently installed in shaft – not removable!

**NOTICE** The appearance and dimensions of the spare parts correspond to the respective nominal volume.

## 12.2.2 Spare parts – volumes of 2.5, 5, and 10 ml

| Vol.            | F + F' | G      | H      | I      | O    |
|-----------------|--------|--------|--------|--------|------|
| 0.25–<br>2.5 ml | 704755 | 704669 | 704689 | 704626 | 7228 |
| 0.5–<br>5 ml    | 704756 | 704606 | 703247 | 704626 | 7288 |
| 1–10 ml         | 704757 | 704607 | 704628 | 704626 | 7298 |

**NOTICE** The appearance and dimensions of the spare parts correspond to the respective nominal volume.

## 12.3 Additional accessories

| Description  | Order No.              |
|--------------|------------------------|
| Label window | <a href="#">704752</a> |
| Label strip  | <a href="#">704753</a> |

## 12 Order Information

| Description  | Order No.              |
|--|------------------------|
| Wall mount   | <a href="#">704882</a> |
| Shelf/rack mount                                     | <a href="#">704881</a> |
| Filter for volume range 2–5 ml, pkg. unit 25 pc      | <a href="#">704652</a> |
| Filter for volume range up to 10 ml, pkg. unit 25 pc | <a href="#">704653</a> |
| Silicone grease for volume range up to 1000 µl       | <a href="#">705502</a> |
| PLT unit (pipette leak detector)                     | <a href="#">703970</a> |

# 13 Repairs

## 13.1 Sending for repair

### NOTICE

Transporting hazardous materials without approval is prohibited by law.

#### Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

#### Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

#### Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

#### Contact addresses

##### Germany:

BRAND GMBH + CO KG  
Otto-Schott-Straße 25  
97877 Wertheim (Germany)

##### USA and Canada:

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)

T +49 9342 808 0  
F +49 9342 808 98000  
[info@brand.de](mailto:info@brand.de)  
[www.brand.de](http://www.brand.de)

T +1-860-767 2562  
F +1 - 860 - 767 2563  
[info@brandtech.com](mailto:info@brandtech.com)  
[www.brandtech.com](http://www.brandtech.com)

**India:**

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park,  
Powai  
Mumbai-400 076 (India)  
T +91 22 42957790  
F +91 22 42957791  
[info@brand.co.in](mailto:info@brand.co.in)  
[www.brand.co.in](http://www.brand.co.in)

**China:**

BRAND (Shanghai) Trading Co., Ltd.  
Rm 201-202, North Tower,  
No. 199 Kaibin Rd, Xuhui District, Shanghai  
Shanghai 200030 (P.R. China)  
T +86 21 6422 2318  
F +86 21 6422 2268  
[info@brand.com.cn](mailto:info@brand.com.cn)  
[www.brand.cn.com](http://www.brand.cn.com)

## 14 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive Liquids.

The complete SOP for testing can be downloaded from [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at [www.brand.de](http://www.brand.de) (Service & Support).

### For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

## 15 Information about your laboratory instrument

The online service MyProduct (<https://www.brand.de/myproduct>) offers quality certificates, accessories, and technical documentation for your laboratory instrument Transferette® pro. By entering the serial and article number there, you will receive information specific to your individual instrument.

On Transferette® pro, you will also find serialized information encoded in the quick response code. This also links to the MyProduct website, where you can access this operating manual and additional certificates for your instrument.

You will also find a data matrix code on some instruments (Transferette® S, HandyStep® touch, and HandyStep touch® S). Scan this using a standard reader app to access the above-mentioned information via the URL <https://www.brand.de/myproduct>.

## 16 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

### **USA and Canada:**

Find more warranty information on [www.brandtech.com](http://www.brandtech.com).

## 17 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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