


## 5 – 100 µl Filter tips 732610 | 732710





Dimensions in mm

Volume	5 – 100 μl
Material	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **

\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.


Status as of: March 2024

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany |T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de



## 5 – 100 µl Filter tips 732630 | 732730





Dimensions in mm

Volume	5 – 100 μl
Material	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP

\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Status as of: March 2024

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany |T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de