


## 50 – 1000 μl Filter tips 732614 | 732714





Dimensions in mm

Volume	50 – 1000 μl
Material	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **

\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml


Status as of: March 2024

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany |T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de



## 50 – 1000 μl Filter tips 732634 | 732734





Dimensions in mm

Volume	50 – 1000 μl
Material	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP
	* human and hastorial DNA _** according to Limulus Amphasute Lysate (LAL) test detection limit 0.01 EU/ml

\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

Status as of: March 2024

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany |T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de