


# Deep well plate 0.3 ml 701355



L	128.15 mm
W	85.70 mm
H	30.20 mm
Hs	2.55 mm
X1	12.10 mm
Xn	115.16 mm
Y1	9.10 mm
Yn	76.60 mm
P	4.50 mm
D	3.65 mm
Hw	26.35 mm
C	24 well
R	16 well

<b>Well Number</b>	384
<b>Type</b>	Standard
<b>Bottom Shape</b>	V-bottom
<b>Nominal Volume</b>	0.3 ml
<b>Material</b>	PP
<b>Well Shape</b>	Square
<b>Cut corner</b>	A1, H1
<b>Quality</b>	BIO-CERT® CERTIFIED QUALITY
<b>Suitable sealing mat</b>	701357



Subject to technical modification without notice. Errors excepted.