



Deep-well plate 1,1 ml Low Profile 701340




L	128.07 mm
W	85.49 mm
H	27.00 mm
Hs	2.22 mm
X1	14.54 mm
Xn	113.53 mm
Y1	11.25 mm
Yn	74.24 mm
P	9.00 mm
C	12 Well
R	8 Well



Well Number	96
Type	Low Profile
Bottom Shape	U-bottom
Nominal Volume	1.1 ml
Material	PP
Well Shape	Round
Cut corner	A12, H12
Quality	CERTIFIED LIFE SCIENCE QUALITY
Suitable sealing mat	701368



Deep-well plate 1,2 ml 701342




L	126.70 mm
W	85.00 mm
H	41.00 mm
Hs	2.95 mm
X1	14.29 mm
Xn	112.41 mm
Y1	11.35 mm
Yn	73.65 mm
P	8.90 mm
C	12 well
R	8 well

Well Number	96
Type	Elevated Rim
Bottom Shape	U-bottom
Nominal Volume	1.2 ml
Material	PP
Well Shape	Round
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY
Suitable sealing mat	701360



Deep-well plate 0,5 ml 701346




L	127.76 mm
W	85.48 mm
H	28.5 mm
Hs	6.10 mm
X1	14.38 mm
Xn	113.38 mm
Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
D	7.20 mm
Hw	21.60 mm
C	12 well
R	8 well

Well Number	96
Type	Standard
Bottom Shape	U-bottom
Nominal Volume	0.5 ml
Material	PP
Well Shape	Square
Cut corner	A1, H1
Quality	BIO-CERT® PCR QUALITY
Suitable sealing mat	701358




Deep-well plate 1 ml U-bottom 701347



Well Number	96
Type	Low Profile
Bottom Shape	U-bottom
Nominal Volume	1 ml
Material	PP
Well Shape	Square
Cut corner	H1
Quality	CERTIFIED LIFE SCIENCE QUALITY
Suitable sealing mat	701362




Deep-well plate 2.2 ml V-bottom, used for King Fisher™ 701348



Well Number	96
Type	Standard
Bottom Shape	V-bottom - pyramid
Nominal Volume	2.2 ml
Material	PP
Well Shape	Square
Cut corner	A1, H1
Quality	CERTIFIED LIFE SCIENCE QUALITY
Suitable sealing mat	701362



Deep-well plate 1.1 ml 701350




L	126.23 mm
W	84.65 mm
H	40.59 mm
Hs	2.97 mm
X1	14.11 mm
Xn	112.12 mm
Y1	11.14 mm
Yn	73.51 mm
P	8.91 mm
D	7.13 mm
Hw	37.03 mm
C	12 well
R	8 well

Well Number	96
Type	Standard
Bottom Shape	U-bottom
Nominal Volume	1.1 ml
Material	PP
Well Shape	Round
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY
Suitable sealing mat	701360



Deep-well plate 2,2 ml V-bottom 701353




L	127.76 mm
W	85.48 mm
H	44.20 mm
Hs	2.60 mm
X1	14.38 mm
Xn	113.38 mm
Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
D	8.24 mm
Hw	39.89 mm
C	12 Well
R	8 Well

Well Number	96
Type	Standard
Bottom Shape	V-bottom
Nominal Volume	2.2 ml
Material	PP
Well Shape	Square
Cut corner	A1, H1
Quality	CERTIFIED LIFE SCIENCE QUALITY
Suitable sealing mat	701362



Deep-well plate 2.2 ml 701354



L	127.76 mm
W	85.47 mm
H	44.00 mm
Hs	3.54 mm
X1	14.32 mm
Xn	113.32 mm
Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
D	8.30 mm
Hw	41.85 mm
C	12 well
R	8 well

Well Number	96
Type	Standard
Bottom Shape	U-bottom
Nominal Volume	2.2 ml
Material	PP
Well Shape	Square
Cut corner	A1, H1
Quality	BIO-CERT® PCR QUALITY
Suitable sealing mat	701362

